



CARE LOCKING IMPLANTS DIMENSIONS



www.orthoimplants.com/ www.orthocare.in
orthocare@vsnl.com/orthocareindia2009@gmail.com

1023

S.No	CAT. NO. Titanium	CAT. NO. Stainless Steel	DESCRIPTION	Thickness	Width	Length
1			LOCKING PLATE CLAVICLE HOOK 3.5mm, RIGHT/LEFT			
	LCP510102	AK-0212R	Locking Plate Clavicle Hook 3.5mm - 4H Right	3.2mm	11.3mm	74mm
	LCP510103	AK-0213R	Locking Plate Clavicle Hook 3.5mm - 5H Right	3.2mm	11.3mm	84mm
	LCP510105	AK-0215R	Locking Plate Clavicle Hook 3.5mm - 6H Right	3.2mm	11.3mm	94mm
	LCP510107	AK-0217R	Locking Plate Clavicle Hook 3.5mm - 7H Right	3.2mm	11.3mm	104mm
	LCP510202	AK-0212L	Locking Plate Clavicle Hook 3.5mm - 4H Left	3.2mm	11.3mm	74mm
	LCP510203	AK-0213L	Locking Plate Clavicle Hook 3.5mm - 5H Left	3.2mm	11.3mm	84mm
	LCP510205	AK-0215L	Locking Plate Clavicle Hook 3.5mm - 6H Left	3.2mm	11.3mm	94mm
	LCP510207	AK-0217L	Locking Plate Clavicle Hook 3.5mm - 7H Left	3.2mm	11.3mm	104mm
2			LOCKING PLATE SUPERIOR ANTERIOR CLAVICLE 2.7/3.5mm WITH LATERAL EXT			
	LCP665303	BC-003R	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 3H Right	3.0mm	9.3mm	69mm
	LCP665304	BC-004R	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 4H Right	3.0mm	9.3mm	81mm
	LCP665305	BC-005R	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 5H Right	3.0mm	9.3mm	94mm
	LCP665306	BC-006R	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 6H Right	3.0mm	9.3mm	105mm
	LCP665307	BC-007R	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 7H Right	3.0mm	9.3mm	123mm
	LCP665308	BC-008R	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 8H Right	3.0mm	9.3mm	135mm
	LCP665403	BC-003L	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 3H Left	3.0mm	9.3mm	69mm
	LCP665404	BC-004L	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 4H Left	3.0mm	9.3mm	81mm
	LCP665405	BC-005L	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 5H Left	3.0mm	9.3mm	94mm
	LCP665406	BC-006L	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 6H Left	3.0mm	9.3mm	105mm
	LCP665407	BC-007L	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 7H Left	3.0mm	9.3mm	123mm
	LCP665408	BC-008L	Locking Plate Superior Anterior Clavicle 2.7/3.5mm With Lateral Extension - 8H Left	3.0mm	9.3mm	135mm
3			LOCKING PLATE SUPERIOR ANTERIOR CLAVICLE 3.5mm, RIGHT/LEFT			
	LCP665106	BA-006R	Locking Plate Superior Anterior Clavicle 3.5mm - 6H Right	3.4mm	10.0mm	94mm
	LCP665107	BA-007R	Locking Plate Superior Anterior Clavicle 3.5mm - 7H Right	3.4mm	10.0mm	110mm
	LCP665108	BA-008R	Locking Plate Superior Anterior Clavicle 3.5mm - 8H Right	3.4mm	10.0mm	120mm
	LCP665206	BA-006L	Locking Plate Superior Anterior Clavicle 3.5mm - 6H Left	3.4mm	10.0mm	94mm

	LCP665207	BA-007L	Locking Plate Superior Anterior Clavicle 3.5mm - 7H Left	3.4mm	10.0mm	110mm
	LCP665208	BA-008L	Locking Plate Superior Anterior Clavicle 3.5mm - 8H Left	3.4mm	10.0mm	120mm
4			LOCKING CLAVICLE S-PLATE 3.5mm RIGHT/LEFT			
	LCPCSP06R	AL-006R	Locking Clavicle S- Plate,6H - RIGHT	2.4mm	10.0mm	85mm
	LCPCSP07R	AL-007R	Locking Clavicle S- Plate,7H - RIGHT	2.4mm	10.0mm	100mm
	LCPCSP08R	AL-008R	Locking Clavicle S- Plate,8H - RIGHT	2.4mm	10.0mm	115mm
	LCPCSP09R	AL-009R	Locking Clavicle S- Plate,9H - RIGHT	2.4mm	10.0mm	130mm
	LCPCSP10R	AL-010R	Locking Clavicle S- Plate,10H - RIGHT	2.4mm	10.0mm	145mm
	LCPCSP12R	AL-012R	Locking Clavicle S- Plate,12H - RIGHT	2.4mm	10.0mm	160mm
	LCPCSP06L	AL-006L	Locking Clavicle S- Plate,6H - LEFT	2.4mm	10.0mm	85mm
	LCPCSP07L	AL-007L	Locking Clavicle S- Plate,7H - LEFT	2.4mm	10.0mm	100mm
	LCPCSP08L	AL-008L	Locking Clavicle S- Plate,8H - LEFT	2.4mm	10.0mm	115mm
	LCPCSP09L	AL-009L	Locking Clavicle S- Plate,9H - LEFT	2.4mm	10.0mm	130mm
	LCPCSP10L	AL-010L	Locking Clavicle S- Plate,10H - LEFT	2.4mm	10.0mm	145mm
	LCPCSP12L	AL-012L	Locking Clavicle S- Plate,12H - LEFT	2.4mm	10.0mm	160mm
5			LOCKING PLATE METAPHYSEAL 3.5mm			
	LCP506006	AD-006	Locking Plate Metaphyseal 3.5mm, 6H	3.2mm	11.0mm	86mm
	LCP506007	AD-007	Locking Plate Metaphyseal 3.5mm, 7H	3.2mm	11.0mm	99mm
	LCP506008	AD-008	Locking Plate Metaphyseal 3.5mm, 8H	3.2mm	11.0mm	112mm
	LCP506009	AD-009	Locking Plate Metaphyseal 3.5mm, 9H	3.2mm	11.0mm	125mm
	LCP506010	AD-010	Locking Plate Metaphyseal 3.5mm, 10H	3.2mm	11.0mm	138mm
	LCP506011	AD-011	Locking Plate Metaphyseal 3.5mm, 11H	3.2mm	11.0mm	151mm
	LCP506012	AD-012	Locking Plate Metaphyseal 3.5mm, 12H	3.2mm	11.0mm	164mm
	LCP506013	AD-013	Locking Plate Metaphyseal 3.5mm, 13H	3.2mm	11.0mm	177mm
	LCP506014	AD-014	Locking Plate Metaphyseal 3.5mm, 14H	3.2mm	11.0mm	190mm
	LCP506015	AD-015	Locking Plate Metaphyseal 3.5mm, 15H	3.2mm	11.0mm	216mm
	LCP506016	AD-016	Locking Plate Metaphyseal 3.5mm, 16H	3.2mm	11.0mm	242mm
	LCP506017	AD-017	Locking Plate Metaphyseal 3.5mm, 17H	3.2mm	11.0mm	260mm
	LCP506018	AD-018	Locking Plate Metaphyseal 3.5mm, 18H	3.2mm	11.0mm	273mm
6			LOCKING PLATE PROXIMAL HUMERAL 3.5mm, SHORT			
	LCP508003	HH-003	Locking Plate Proximal Humeral 3.5mm, Short - 3H	3.2mm	13.0mm	90mm
	LCP508005	HH-005	Locking Plate Proximal Humeral 3.5mm, Short - 5H	3.2mm	13.0mm	115mm
7			LOCKING PLATE PROXIMAL HUMERAL 3.5mm, LONG			
	LCP509205	GG-005	Locking Plate Proximal Humeral 3.5mm, Long - 5H	3.2mm	13.0mm	142mm
	LCP509206	GG-006	Locking Plate Proximal Humeral 3.5mm, Long - 6H	3.2mm	13.0mm	160mm
	LCP509208	GG-008	Locking Plate Proximal Humeral 3.5mm, Long - 8H	3.2mm	13.0mm	196mm
	LCP509209	GG-009	Locking Plate Proximal Humeral 3.5mm, Long - 9H	3.2mm	13.0mm	224mm

	LCP509210	GG-010	Locking Plate Proximal Humeral 3.5mm, Long - 10H	3.2mm	13.0mm	232mm
	LCP509212	GG-012	Locking Plate Proximal Humeral 3.5mm, Long - 12H	3.2mm	13.0mm	268mm
8			LOCKING PLATE PROXIMAL HUMERAL 3.5mm TYPE - I			
	LCP519205	PH-005	Locking Plate Proximal Humeral 3.5mm , Type I -5H	3.2mm	13.0mm	102mm
	LCP519206	PH-006	Locking Plate Proximal Humeral 3.5mm , Type I -6H	3.2mm	13.0mm	120mm
	LCP519208	PH-008	Locking Plate Proximal Humeral 3.5mm , Type I -8H	3.2mm	13.0mm	138mm
	LCP519209	PH-009	Locking Plate Proximal Humeral 3.5mm , Type I -9H	3.2mm	13.0mm	156mm
	LCP519210	PH-010	Locking Plate Proximal Humeral 3.5mm , Type I -10H	3.2mm	13.0mm	174mm
	LCP519212	PH-012	Locking Plate Proximal Humeral 3.5mm , Type I -12H	3.2mm	13.0mm	192mm
9			LOCKING PLATE PERIARTICULAR PROXIMAL HUMERAL 3.5mm RIGHT/LEFT			
	LCP681102	PPHP-002R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 2H, Right	3.2mm	13.0mm	91mm
	LCP681103	PPHP-003R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 3H, Right	3.2mm	13.0mm	127mm
	LCP681104	PPHP-004R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 4H, Right	3.2mm	13.0mm	163mm
	LCP681105	PPHP-005R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 5H, Right	3.2mm	13.0mm	199mm
	LCP681106	PPHP-006R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 6H, Right	3.2mm	13.0mm	235mm
	LCP681108	PPHP-008R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 8H, Right	3.2mm	13.0mm	271mm
	LCP681110	PPHP-010R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 10H, Right	3.2mm	13.0mm	307mm
	LCP681112	PPHP-012R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 12H, Right	3.2mm	13.0mm	343mm
	LCP681114	PPHP-014R	Locking Plate Periarticular Proximal Humeral 3.5 mm, 14H, Right	3.2mm	13.0mm	379mm
	LCP681202	PPHP-002L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 2H, Left	3.2mm	13.0mm	91mm
	LCP681203	PPHP-003L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 3H, Left	3.2mm	13.0mm	127mm
	LCP681204	PPHP-004L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 4H, Left	3.2mm	13.0mm	163mm
	LCP681205	PPHP-005L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 5H, Left	3.2mm	13.0mm	199mm
	LCP681206	PPHP-006L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 6H, Left	3.2mm	13.0mm	235mm
	LCP681208	PPHP-008L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 8H, Left	3.2mm	13.0mm	271mm
	LCP681210	PPHP-010L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 10H, Left	3.2mm	13.0mm	307mm
	LCP681212	PPHP-012L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 12H, Left	3.2mm	13.0mm	343mm
	LCP681214	PPHP-014L	Locking Plate Periarticular Proximal Humeral 3.5 mm, 14H, Left	3.2mm	13.0mm	379mm
10			LOCKING PLATE CLOVER LEAF 3.5mm			
	LCP676903	BD-003	Locking Plate Clover Leaf 3.5mm - 3H	8.0mm	15.8mm	88mm
	LCP676904	BD-004	Locking Plate Clover Leaf 3.5mm - 4H	8.0mm	15.8mm	104mm
	LCP676905	BD-005	Locking Plate Clover Leaf 3.5mm - 5H	8.0mm	15.8mm	120mm
	LCP676906	BD-006	Locking Plate Clover Leaf 3.5mm - 6H	8.0mm	15.8mm	136mm
	LCP676907	BD-007	Locking Plate Clover Leaf 3.5mm - 7H	8.0mm	15.8mm	152mm
	LCP676908	BD-008	Locking Plate Clover Leaf 3.5mm - 8H	8.0mm	15.8mm	168mm
11			LOCKING PLATE SMALL 3.5mm			
	LCP501005	AA-005	Locking Plate Small 3.5mm - 5H	3.3mm	11.0mm	72mm
	LCP501006	AA-006	Locking Plate Small 3.5mm - 6H	3.3mm	11.0mm	85mm

	LCP501007	AA-007	Locking Plate Small 3.5mm - 7H	3.3mm	11.0mm	98mm
	LCP501008	AA-008	Locking Plate Small 3.5mm - 8H	3.3mm	11.0mm	111mm
	LCP501009	AA-009	Locking Plate Small 3.5mm - 9H	3.3mm	11.0mm	124mm
	LCP501010	AA-010	Locking Plate Small 3.5mm - 10H	3.3mm	11.0mm	137mm
	LCP501011	AA-011	Locking Plate Small 3.5mm - 11H	3.3mm	11.0mm	150mm
	LCP501012	AA-012	Locking Plate Small 3.5mm - 12H	3.3mm	11.0mm	163mm
12			LOCKING PLATE CURVED NARROW 3.5mm			
	LCP700010	CN-010	Locking Plate Curved Narrow 3.5mm - 10H	4.5mm	13.5mm	150mm
	LCP700011	CN-011	Locking Plate Curved Narrow 3.5mm - 11H	4.5mm	13.5mm	163mm
	LCP700012	CN-012	Locking Plate Curved Narrow 3.5mm - 12H	4.5mm	13.5mm	176mm
	LCP700014	CN-014	Locking Plate Curved Narrow 3.5mm - 14H	4.5mm	13.5mm	189mm
	LCP700016	CN-016	Locking Plate Curved Narrow 3.5mm - 16H	4.5mm	13.5mm	202mm
	LCP700018	CN-018	Locking Plate Curved Narrow 3.5mm - 18H	4.5mm	13.5mm	215mm
	LCP700020	CN-020	Locking Plate Curved Narrow 3.5mm - 20H	4.5mm	13.5mm	228mm
	LCP700022	CN-022	Locking Plate Curved Narrow 3.5mm - 22H	4.5mm	13.5mm	241mm
13			LOCKING PLATE CURVED BROAD 3.5mm			
	LCP710010	CB-010	Locking Plate Curved Broad 3.5mm - 10H	5.6mm	17.5mm	150mm
	LCP710011	CB-011	Locking Plate Curved Broad 3.5mm - 11H	5.6mm	17.5mm	163mm
	LCP710012	CB-012	Locking Plate Curved Broad 3.5mm - 12H	5.6mm	17.5mm	176mm
	LCP710014	CB-014	Locking Plate Curved Broad 3.5mm - 14H	5.6mm	17.5mm	189mm
	LCP710016	CB-016	Locking Plate Curved Broad 3.5mm - 16H	5.6mm	17.5mm	202mm
	LCP710018	CB-018	Locking Plate Curved Broad 3.5mm - 18H	5.6mm	17.5mm	215mm
	LCP710020	CB-020	Locking Plate Curved Broad 3.5mm - 20H	5.6mm	17.5mm	228mm
	LCP710022	CB-022	Locking Plate Curved Broad 3.5mm - 22H	5.6mm	17.5mm	241mm
	LCP710024	CB-024	Locking Plate Curved Broad 3.5mm - 24H	5.6mm	17.5mm	254mm
	LCP710026	CB-026	Locking Plate Curved Broad 3.5mm - 26H	5.6mm	17.5mm	267mm
	LCP710028	CB-028	Locking Plate Curved Broad 3.5mm - 28H	5.6mm	17.5mm	280mm
	LCP710030	CB-030	Locking Plate Curved Broad 3.5mm - 30H	5.6mm	17.5mm	293mm
14	LCP666503	BJ-003	LOCKING HOOK PLATE 3.5mm, 3 Holes Length - 62mm	1.0mm	4.0mm	62mm
15			LOCKING PLATE ONE -THIRD TUBULAR 3.5mm			
	LCP681002	OT-002	Locking One-Third Tubular Plate 3.5mm - 2 Holes	1.8mm	10.0mm	28mm
	LCP681003	OT-003	Locking One-Third Tubular Plate 3.5mm - 3 Holes	1.8mm	10.0mm	40mm
	LCP681004	OT-004	Locking One-Third Tubular Plate 3.5mm - 4 Holes	1.8mm	10.0mm	52mm
	LCP681005	OT-005	Locking One-Third Tubular Plate 3.5mm - 5 Holes	1.8mm	10.0mm	64mm
	LCP681006	OT-006	Locking One-Third Tubular Plate 3.5mm - 6 Holes	1.8mm	10.0mm	76mm
	LCP681007	OT-007	Locking One-Third Tubular Plate 3.5mm - 7 Holes	1.8mm	10.0mm	88mm
	LCP681008	OT-008	Locking One-Third Tubular Plate 3.5mm - 8 Holes	1.8mm	10.0mm	100mm
	LCP681009	OT-009	Locking One-Third Tubular Plate 3.5mm - 9 Holes	1.8mm	10.0mm	112mm
	LCP681010	OT-010	Locking One-Third Tubular Plate 3.5mm - 10 Holes	1.8mm	10.0mm	124mm

	LCP681012	OT-012	Locking One-Third Tubular Plate 3.5mm - 12 Holes	1.8mm	10.0mm	148mm
16			LOCKING PLATE HUMERAL DISTAL MEDIAL 2.7/3.5mm, RIGHT/LEFT			
	LCP513203	AY-003R	Locking Plate Humeral Distal Medial 2.7/3.5mm - 3H Right	3.8mm	11.0mm	59mm
	LCP513205	AY-005R	Locking Plate Humeral Distal Medial 2.7/3.5mm - 5H Right	3.8mm	11.0mm	84mm
	LCP513207	AY-007R	Locking Plate Humeral Distal Medial 2.7/3.5mm - 7H Right	3.8mm	11.0mm	110mm
	LCP513209	AY-009R	Locking Plate Humeral Distal Medial 2.7/3.5mm - 9H Right	3.8mm	11.0mm	136mm
	LCP513211	AY-011R	Locking Plate Humeral Distal Medial 2.7/3.5mm - 11H Right	3.8mm	11.0mm	201mm
	LCP513103	AY-003L	Locking Plate Humeral Distal Medial 2.7/3.5mm - 3H Left	3.8mm	11.0mm	59mm
	LCP513105	AY-005L	Locking Plate Humeral Distal Medial 2.7/3.5mm - 5H Left	3.8mm	11.0mm	84mm
	LCP513107	AY-007L	Locking Plate Humeral Distal Medial 2.7/3.5mm - 7H Left	3.8mm	11.0mm	110mm
	LCP513109	AY-009L	Locking Plate Humeral Distal Medial 2.7/3.5mm - 9H Left	3.8mm	11.0mm	136mm
	LCP513111	AY-011L	Locking Plate Humeral Distal Medial 2.7/3.5mm - 11H Left	3.8mm	11.0mm	201mm
17			LOCKING PLATE HUMERAL DISTAL MEDIAL 3.5mm			
	LCP513403	AI-003	Locking Plate Humeral Distal Medial 3.5mm - 3H	3.8mm	11.0mm	59mm
	LCP513405	AI-005	Locking Plate Humeral Distal Medial 3.5mm - 5H	3.8mm	11.0mm	84mm
	LCP513407	AI-007	Locking Plate Humeral Distal Medial 3.5mm - 7H	3.8mm	11.0mm	110mm
	LCP513409	AI-009	Locking Plate Humeral Distal Medial 3.5mm - 9H	3.8mm	11.0mm	136mm
	LCP513411	AI-011	Locking Plate Humeral Distal Medial 3.5mm - 11H	3.8mm	11.0mm	201mm
18			LOCKING PLATE HUMERAL DISTAL DORSAL WITHOUT SUPPORT 2.7/3.5mm, R/L			
	LCP504403	AZ-003R	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 3H Right	3.8mm	11.0mm	65mm
	LCP504405	AZ-005R	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 5H Right	3.8mm	11.0mm	91mm
	LCP504407	AZ-007R	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 7H Right	3.8mm	11.0mm	117mm
	LCP504409	AZ-009R	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 9H Right	3.8mm	11.0mm	143mm
	LCP504411	AZ-011R	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 11H Right	3.8mm	11.0mm	169mm
	LCP504414	AZ-014R	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 14H Right	3.8mm	11.0mm	208mm
	LCP504303	AZ-003L	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 3H Left	3.8mm	11.0mm	65mm
	LCP504305	AZ-005L	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 5H Left	3.8mm	11.0mm	91mm
	LCP504307	AZ-007L	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 7H Left	3.8mm	11.0mm	117mm
	LCP504309	AZ-009L	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 9H Left	3.8mm	11.0mm	143mm
	LCP504311	AZ-011L	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 11H Left	3.8mm	11.0mm	169mm
	LCP504314	AZ-014L	Locking Plate Humeral Distal Dorsal Without Support 2.7/3.5mm - 14H Left	3.8mm	11.0mm	208mm
19			LOCKING PLATE HUMERAL DISTAL DORSAL WITHOUT SUPPORT 3.5mm, R/L			
	LCP504703	AM-003R	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 3H Right	3.8mm	11.0mm	65mm
	LCP504705	AM-005R	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 5H Right	3.8mm	11.0mm	91mm
	LCP504707	AM-007R	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 7H Right	3.8mm	11.0mm	117mm
	LCP504709	AM-009R	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 9H Right	3.8mm	11.0mm	143mm
	LCP504711	AM-011R	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 11H Right	3.8mm	11.0mm	169mm

	LCP504714	AM-014R	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 14H Right	3.8mm	11.0mm	208mm	
	LCP504803	AM-003L	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 3H Left	3.8mm	11.0mm	65mm	
	LCP504805	AM-005L	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 5H Left	3.8mm	11.0mm	91mm	
	LCP504807	AM-007L	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 7H Left	3.8mm	11.0mm	117mm	
	LCP504809	AM-009L	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 9H Left	3.8mm	11.0mm	143mm	
	LCP504811	AM-011L	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 11H Left	3.8mm	11.0mm	169mm	
	LCP504814	AM-014L	Locking Plate Humeral Distal Dorsal Without Support 3.5mm - 14H Left	3.8mm	11.0mm	208mm	
20			LOCKING PLATE HUMERAL DISTAL DORSAL WITH SUPPORT 2.7/3.5mm, R/L				
	LCP504503	AK-003R	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 3H Right	3.8mm	11.0mm	65mm	
	LCP504505	AK-005R	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 5H Right	3.8mm	11.0mm	91mm	
	LCP504507	AK-007R	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 7H Right	3.8mm	11.0mm	117mm	
	LCP504509	AK-009R	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 9H Right	3.8mm	11.0mm	143mm	
	LCP504511	AK-011R	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 11H Right	3.8mm	11.0mm	169mm	
	LCP504514	AK-014R	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 14H Right	3.8mm	11.0mm	208mm	
	LCP504603	AK-003L	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 3H Left	3.8mm	11.0mm	65mm	
	LCP504605	AK-005L	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 5H Left	3.8mm	11.0mm	91mm	
	LCP504607	AK-007L	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 7H Left	3.8mm	11.0mm	117mm	
	LCP504609	AK-009L	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 9H Left	3.8mm	11.0mm	143mm	
	LCP504611	AK-011L	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 11H Left	3.8mm	11.0mm	169mm	
	LCP504614	AK-014L	Locking Plate Humeral Distal Dorsal With Support 2.7/3.5mm - 14H Left	3.8mm	11.0mm	208mm	
21			LOCKING PLATE HUMERAL DISTAL DORSAL WITH SUPPORT 3.5mm RIGHT/LEFT				
	LCP504103	NOT AVAILABLE	Locking Plate Humeral Distal Dorsal With Support 3.5mm - 3H Right	3.8mm	11.0mm	65mm	
	LCP504105		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 5H Right	3.8mm	11.0mm	91mm	
	LCP504107		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 7H Right	3.8mm	11.0mm	117mm	
	LCP504109		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 9H Right	3.8mm	11.0mm	143mm	
	LCP504211		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 11H Right	3.8mm	11.0mm	169mm	
	LCP504214		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 14H Right	3.8mm	11.0mm	208mm	
	LCP504203		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 3H Left	3.8mm	11.0mm	65mm	
	LCP504205		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 5H Left	3.8mm	11.0mm	91mm	
	LCP504207		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 7H Left	3.8mm	11.0mm	117mm	
	LCP504209		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 9H Left	3.8mm	11.0mm	143mm	
	LCP504211		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 11H Left	3.8mm	11.0mm	169mm	
	LCP504214		Locking Plate Humeral Distal Dorsal With Support 3.5mm - 14H Left	3.8mm	11.0mm	208mm	
22				LOCKING PLATE EXTRA-ARTICULAR DISTAL HUMERUS 3.5mm, RIGHT/LEFT			
	LCP666204		BH-004R	Locking Plate Extra-Articular Distal Humerus 3.5mm - 4H Right	4.3mm	12.7mm	122mm
	LCP666206	BH-006R	Locking Plate Extra-Articular Distal Humerus 3.5mm - 6H Right	4.3mm	12.7mm	158mm	
	LCP666208	BH-008R	Locking Plate Extra-Articular Distal Humerus 3.5mm - 8H Right	4.3mm	12.7mm	194mm	

	LCP666210	BH-010R	Locking Plate Extra-Articular Distal Humerus 3.5mm - 10H Right	4.3mm	12.7mm	230mm
	LCP666212	BH-012R	Locking Plate Extra-Articular Distal Humerus 3.5mm - 12H Right	4.3mm	12.7mm	266mm
	LCP666214	BH-014R	Locking Plate Extra-Articular Distal Humerus 3.5mm - 14H Right	4.3mm	12.7mm	302mm
	LCP666104	BH-004L	Locking Plate Extra-Articular Distal Humerus 3.5mm - 4H Left	4.3mm	12.7mm	122mm
	LCP666106	BH-006L	Locking Plate Extra-Articular Distal Humerus 3.5mm - 6H Left	4.3mm	12.7mm	158mm
	LCP666108	BH-008L	Locking Plate Extra-Articular Distal Humerus 3.5mm - 8H Left	4.3mm	12.7mm	194mm
	LCP666110	BH-010L	Locking Plate Extra-Articular Distal Humerus 3.5mm - 10H Left	4.3mm	12.7mm	230mm
	LCP666112	BH-012L	Locking Plate Extra-Articular Distal Humerus 3.5mm - 12H Left	4.3mm	12.7mm	266mm
	LCP666114	BH-014L	Locking Plate Extra-Articular Distal Humerus 3.5mm - 14H Left	4.3mm	12.7mm	302mm
23			LOCKING PLATE OLECRANON 3.5mm TYPE - I, RIGHT/LEFT			
	LCP666302	BI-002R	Locking Plate Olecranon 3.5mm Type-I - 2H Right	2.4mm	10.0mm	86mm
	LCP666304	BI-004R	Locking Plate Olecranon 3.5mm Type-I - 4H Right	2.4mm	10.0mm	112mm
	LCP666306	BI-006R	Locking Plate Olecranon 3.5mm Type-I - 6H Right	2.4mm	10.0mm	137mm
	LCP666308	BI-008R	Locking Plate Olecranon 3.5mm Type-I - 8H Right	2.4mm	10.0mm	163mm
	LCP666310	BI-010R	Locking Plate Olecranon 3.5mm Type-I - 10H Right	2.4mm	10.0mm	189mm
	LCP666312	BI-012R	Locking Plate Olecranon 3.5mm Type-I - 12H Right	2.4mm	10.0mm	215mm
	LCP666402	BI-002L	Locking Plate Olecranon 3.5mm Type-I - 2H Left	2.4mm	10.0mm	86mm
	LCP666404	BI-004L	Locking Plate Olecranon 3.5mm Type-I - 4H Left	2.4mm	10.0mm	112mm
	LCP666406	BI-006L	Locking Plate Olecranon 3.5mm Type-I - 6H Left	2.4mm	10.0mm	137mm
	LCP666408	BI-008L	Locking Plate Olecranon 3.5mm Type-I - 8H Left	2.4mm	10.0mm	163mm
	LCP666410	BI-010L	Locking Plate Olecranon 3.5mm Type-I - 10H Left	2.4mm	10.0mm	189mm
	LCP666412	BI-012L	Locking Plate Olecranon 3.5mm Type-I - 12H Left	2.4mm	10.0mm	215mm
24			LOCKING PLATE OLECRANON 3.5mm TYPE - II			
	LCP515003	AH-003	Locking Plate Olecranon 3.5mm Type-II - 3H	2.8mm	10.0mm	88mm
	LCP515004	AH-004	Locking Plate Olecranon 3.5mm Type-II - 4H	2.8mm	10.0mm	96mm
	LCP515005	AH-005	Locking Plate Olecranon 3.5mm Type-II - 5H	2.8mm	10.0mm	104mm
	LCP515006	AH-006	Locking Plate Olecranon 3.5mm Type-II - 6H	2.8mm	10.0mm	112mm
	LCP515007	AH-007	Locking Plate Olecranon 3.5mm Type-II - 7H	2.8mm	10.0mm	120mm
	LCP515008	AH-008	Locking Plate Olecranon 3.5mm Type-II - 8H	2.8mm	10.0mm	128mm
25			LOCKING T-PLATE RIGHT ANGLE 3.5mm			
	LCP512003	BB-003	Locking T-Plate Right Angle 3.5mm - 3H	2.0mm	10.0mm	50mm
	LCP512004	BB-004	Locking T-Plate Right Angle 3.5mm - 4H	2.0mm	10.0mm	59mm
	LCP512005	BB-005	Locking T-Plate Right Angle 3.5mm - 5H	2.0mm	10.0mm	67mm
	LCP512006	BB-006	Locking T-Plate Right Angle 3.5mm - 6H	2.0mm	10.0mm	76mm
	LCP512007	BB-007	Locking T-Plate Right Angle 3.5mm - 7H	2.0mm	10.0mm	85mm
26			LOCKING T-PLATE RIGHT ANGLE 3.5mm, 4 HOLE HEAD			

	LCP514003	LCBB003	Locking T-Plate Right Angle 3.5mm, 4 Hole Head - 3H	2.0mm	10.0mm	50mm
	LCP514004	LCBB004	Locking T-Plate Right Angle 3.5mm, 4 Hole Head - 4H	2.0mm	10.0mm	59mm
	LCP514005	LCBB005	Locking T-Plate Right Angle 3.5mm, 4 Hole Head - 5H	2.0mm	10.0mm	67mm
	LCP514006	LCBB006	Locking T-Plate Right Angle 3.5mm, 4 Hole Head - 6H	2.0mm	10.0mm	76mm
				2.0mm	10.0mm	85mm
27			LOCKING T-PLATE OBLIQUE ANGLE 3.5mm, RIGHT/LEFT			
	LCP511103	CC-003R	Locking T-Plate Oblique Angle 3.5mm - 3H Right	2.0mm	10.0mm	52mm
	LCP511104	CC-004R	Locking T-Plate Oblique Angle 3.5mm - 4H Right	2.0mm	10.0mm	63mm
	LCP511105	CC-005R	Locking T-Plate Oblique Angle 3.5mm - 5H Right	2.0mm	10.0mm	74mm
	LCP511203	CC-003L	Locking T-Plate Oblique Angle 3.5mm - 3H Left	2.0mm	10.0mm	52mm
	LCP511204	CC-004L	Locking T-Plate Oblique Angle 3.5mm - 4H Left	2.0mm	10.0mm	63mm
	LCP511205	CC-005L	Locking T-Plate Oblique Angle 3.5mm - 5H Left	2.0mm	10.0mm	74mm
28			LOCKING DISTAL VOLAR RADIUS ANATOMIC PLATES 2.4/2.7mm , R/L			
	LCP677003	DVRS-003R	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 3H Right	2.25mm	8.2mm	49mm
	LCP677004	DVRS-004R	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 4H Right	2.25mm	8.2mm	58mm
	LCP677007	DVRS-007R	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 7H Right	2.25mm	8.2mm	85mm
	LCP677012	DVRS-012R	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 12H Right	2.25mm	8.2mm	130mm
	LCP677103	DVRS-003L	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 3H Left	2.25mm	8.2mm	49mm
	LCP677104	DVRS-004L	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 4H Left	2.25mm	8.2mm	58mm
	LCP677107	DVRS-007L	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 7H Left	2.25mm	8.2mm	85mm
	LCP677112	DVRS-012L	Locking Distal Volar Radius Anatomic Plates 2.4/2.7mm, - 12H Left	2.25mm	8.2mm	130mm
29			LOCKING PLATE VOLAR COLUMN DISTAL RADIUS 2.4/2.7mm			
			LOCKING PLATE VOLAR COLUMN DISTAL RADIUS 2.4/2.7mm 8 HOLE HEAD, R/L			
	LCP666903	BL-003R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 8 Hole Head - 3H Right	2.25mm	8.2mm	49mm
	LCP666904	BL-004R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 8 Hole Head - 4H Right	2.25mm	8.2mm	58mm
	LCP666905	BL-005R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 8 Hole Head - 5H Right	2.25mm	8.2mm	67mm
	LCP666803	BL-003L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 8 Hole Head - 3H Left	2.25mm	8.2mm	49mm
	LCP666804	BL-004L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 8 Hole Head - 4H Left	2.25mm	8.2mm	58mm
	LCP666805	BL-005L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 8 Hole Head - 5H Left	2.25mm	8.2mm	67mm
			LOCKING PLATE VOLAR COLUMN DISTAL RADIUS 2.4/2.7mm 9 HOLE HEAD, R/L			
	LCP666703	BK-003R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 9 Hole Head - 3H Right	2.25mm	8.2mm	49mm
	LCP666704	BK-004R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 9 Hole Head - 4H Right	2.25mm	8.2mm	58mm
	LCP666705	BK-005R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 9 Hole Head - 5H Right	2.25mm	8.2mm	67mm
	LCP666603	BK-003L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 9 Hole Head - 3H Left	2.25mm	8.2mm	49mm
	LCP666604	BK-004L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 9 Hole Head - 4H Left	2.25mm	8.2mm	58mm
	LCP666605	BK-005L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 9 Hole Head - 5H Left	2.25mm	8.2mm	67mm

30			LOCKING PLATE PROXIMAL RADIUS 2.4/2.7mm			
			LOCKING RIM PLATE FOR RADIAL HEAD RIGHT/LEFT			
	LCP720102	RHP-002R	Locking Rim Plate for Radial Head - 2H Right	1.8mm	8.0mm	36mm
	LCP720103	RHP-003R	Locking Rim Plate for Radial Head - 3H Right	1.8mm	8.0mm	46mm
	LCP720104	RHP-004R	Locking Rim Plate for Radial Head - 4H Right	1.8mm	8.0mm	56mm
	LCP720202	RHP-002L	Locking Rim Plate for Radial Head - 2H Left	1.8mm	8.0mm	36mm
	LCP720203	RHP-003L	Locking Rim Plate for Radial Head - 3H Left	1.8mm	8.0mm	46mm
	LCP720204	RHP-004L	Locking Rim Plate for Radial Head - 4H Left	1.8mm	8.0mm	56mm
			LOCKING NECK PLATE FOR RADIAL HEAD			
	LCP730002	NHP-002	Locking Neck Plate for Radial Head - 2H	1.8mm	8.0mm	30mm
	LCP730003	NHP-003	Locking Neck Plate for Radial Head - 3H	1.8mm	8.0mm	40mm
	LCP730004	NHP-004	Locking Neck Plate for Radial Head - 4H	1.8mm	8.0mm	50mm
31			LOCKING PLATE VOLAR COLUMN DISTAL RADIUS 2.4/2.7mm 6 HOLE HEAD, R/L			
			LOCKING PLATE VOLAR COLUMN DISTAL RADIUS 2.4/2.7mm VARIABLE ANGLE, R/L			
	LCP691602	VC-002R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 2H Right	2.25mm	8.2mm	36mm
	LCP691603	VC-003R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 3H Right	2.25mm	8.2mm	45mm
	LCP691604	VC-004R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 4H Right	2.25mm	8.2mm	54mm
	LCP691605	VC-005R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 5H Right	2.25mm	8.2mm	63mm
	LCP691502	VC-002L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 2H Left	2.25mm	8.2mm	36mm
	LCP691503	VC-003L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 3H Left	2.25mm	8.2mm	45mm
	LCP691504	VC-004L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 4H Left	2.25mm	8.2mm	54mm
	LCP691505	VC-005L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, 6 Hole Head - 5H Left	2.25mm	8.2mm	63mm
			LOCKING PLATE VOLAR COLUMN DISTAL RADIUS 2.4/2.7mm FIX ANGLE R/L			
	LCPV691602	VAC-002R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 2H Right	2.25mm	8.2mm	36mm
	LCPV691603	VAC-003R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 3H Right	2.25mm	8.2mm	45mm
	LCPV691604	VAC-004R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 4H Right	2.25mm	8.2mm	54mm
	LCPV691605	VAC-005R	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 5H Right	2.25mm	8.2mm	63mm
	LCPV691502	VAC-002L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 2H Left	2.25mm	8.2mm	36mm
	LCPV691503	VAC-003L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 3H Left	2.25mm	8.2mm	45mm
	LCPV691504	VAC-004L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 4H Left	2.25mm	8.2mm	54mm
	LCPV691505	VAC-005L	Locking Plate Volar Column Distal Radius 2.4/2.7mm, Fix Angel- 5H Left	2.25mm	8.2mm	63mm
32			LOCKING PLATE EXTRA-ARTICULAR VOLAR DISTAL RADIUS 2.4/2.7mm 4 HOLE H			
	LCP676103	VA-003R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 4 Hole Head - 3H Right	2.6mm	7.0mm	47mm
	LCP676104	VA-004R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 4 Hole Head - 4H Right	2.6mm	7.0mm	56mm
	LCP676105	VA-005R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 4 Hole Head - 5H Right	2.6mm	7.0mm	65mm

	LCP676203	VA-003L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 4 Hole Head - 3H Left	2.6mm	7.0mm	47mm
	LCP676204	VA-004L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 4 Hole Head - 4H Left	2.6mm	7.0mm	56mm
	LCP676205	VA-005L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 4 Hole Head - 5H Left	2.6mm	7.0mm	65mm
33			LOCKING PLATE EXTRA-ARTICULAR VOLAR DISTAL RADIUS 2.4/2.7mm 5 HOLE H			
	LCP676303	VD-003R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- I - 3H Right	2.6mm	7.0mm	47mm
	LCP676304	VD-004R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- I - 4H Right	2.6mm	7.0mm	56mm
	LCP676305	VD-005R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- I - 5H Right	2.6mm	7.0mm	65mm
	LCP676403	VD-003L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- I - 3H Left	2.6mm	7.0mm	47mm
	LCP676404	VD-004L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- I - 4H Left	2.6mm	7.0mm	56mm
	LCP676405	VD-005L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- I - 5H Left	2.6mm	7.0mm	65mm
34			LOCKING PLATE EXTRA-ARTICULAR VOLAR DISTAL RADIUS 2.4/2.7mm 5 HOLE H			
	LCP676503	VR-003R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- II - 3H Right	2.6mm	7.0mm	43mm
	LCP676504	VR-004R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- II - 4H Right	2.6mm	7.0mm	52mm
	LCP676505	VR-005R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- II - 5H Right	2.6mm	7.0mm	61mm
	LCP676603	VR-003L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- II - 3H Left	2.6mm	7.0mm	43mm
	LCP676604	VR-004L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- II - 4H Left	2.6mm	7.0mm	52mm
	LCP676605	VR-005L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- II - 5H Left	2.6mm	7.0mm	61mm
35			LOCKING PLATE EXTRA-ARTICULAR VOLAR DISTAL RADIUS 2.4/2.7mm 5 HOLE H			
	LCP676703	VS-003R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type-III - 3H Right	2.6mm	7.0mm	47mm
	LCP676704	VS-004R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- III - 4H Right	2.6mm	7.0mm	56mm
	LCP676705	VS-005R	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- III - 5H Right	2.6mm	7.0mm	65mm
	LCP676803	VS-003L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- III - 3H Left	2.6mm	7.0mm	47mm
	LCP676804	VS-004L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- III - 4H Left	2.6mm	7.0mm	56mm
	LCP676805	VS-005L	Locking Plate Extra-Articular Volar Distal Radius 2.4/2.7mm, 5 Hole Head Type- III - 5H Left	2.6mm	7.0mm	65mm
36			LOCKING PLATE WRIST FUSION			
	LCP667001	BM-000	Locking Plate Wrist Fusion Standard Bend	3.4mm	11.0mm	118mm
	LCP667002	BN-000	Locking Plate Wrist Fusion Short Bend	3.4mm	11.0mm	118mm
	LCP667003	BO-000	Locking Plate Wrist Fusion Straight Bend	3.4mm	11.0mm	118mm
37			LOCKING PLATE DISTAL LATERAL RADIUS 2.7mm			
	LCP556005	DLR-6005	Locking Plate Distal Lateral Radius 2.7mm - 5H	1.6mm	6.5mm	46mm
	LCP556006	DLR-6006	Locking Plate Distal Lateral Radius 2.7mm - 6H	1.6mm	6.5mm	57mm
38			LOCKING DISTAL DORSAL RADIUS T-PLATE 2.7mm			
	LCP557003	DDR-7003	Locking Distal Dorsal Radius T-Plate 2.7mm - 3H	1.6mm	6.4mm	37mm

	LCP557004	DDR-7004	Locking Distal Dorsal Radius T-Plate 2.7mm - 4H	1.6mm	6.4mm	44mm
39			LOCKING DISTAL DORSAL RADIUS L-PLATE 2.7mm, RIGHT/LEFT			
	LCP556103	DDR-6103	Locking Distal Dorsal Radius L-Plate 2.7mm - 3H Right	1.6mm	6.4mm	37mm
	LCP556104	DDR-6104	Locking Distal Dorsal Radius L-Plate 2.7mm - 4H Right	1.6mm	6.4mm	44mm
	LCP556203	DDR-6203	Locking Distal Dorsal Radius L-Plate 2.7mm - 3H Left	1.6mm	6.4mm	37mm
	LCP556204	DDR-6204	Locking Distal Dorsal Radius L-Plate 2.7mm - 4H Left	1.6mm	6.4mm	44mm
40			LOCKING DISTAL DORSAL RADIUS L-PLATE 2.7mm Type-II, RIGHT/LEFT			
	LCP558103	DDR-8103	Locking Distal Dorsal Radius L-Plate 2.7mm, Type-II - 3H Right	1.6mm	6.4mm	41mm
	LCP558104	DDR-8104	Locking Distal Dorsal Radius L-Plate 2.7mm, Type-II - 4H Right	1.6mm	6.4mm	49mm
	LCP558203	DDR-8203	Locking Distal Dorsal Radius L-Plate 2.7mm, Type-II - 3H Left	1.6mm	6.4mm	41mm
	LCP558204	DDR-8204	Locking Distal Dorsal Radius L-Plate 2.7mm, Type-II - 4H Left	1.6mm	6.4mm	49mm
41			LOCKING PLATE RADIO DISTAL DORSAL 3.5mm, RIGHT/LEFT			
	LCP503103	AG-003R	Locking Plate Radio Distal Dorsal 3.5mm - 3H Right	2.0mm	10.0mm	48mm
	LCP503104	AG-004R	Locking Plate Radio Distal Dorsal 3.5mm - 4H Right	2.0mm	10.0mm	60mm
	LCP503106	AG-006R	Locking Plate Radio Distal Dorsal 3.5mm - 6H Right	2.0mm	10.0mm	84mm
	LCP503108	AG-008R	Locking Plate Radio Distal Dorsal 3.5mm - 8H Right	2.0mm	10.0mm	108mm
	LCP503110	AG-010R	Locking Plate Radio Distal Dorsal 3.5mm - 10H Right	2.0mm	10.0mm	132mm
	LCP503203	AG-003L	Locking Plate Radio Distal Dorsal 3.5mm - 3H Left	2.0mm	10.0mm	48mm
	LCP503204	AG-004L	Locking Plate Radio Distal Dorsal 3.5mm - 4H Left	2.0mm	10.0mm	60mm
	LCP503206	AG-006L	Locking Plate Radio Distal Dorsal 3.5mm - 6H Left	2.0mm	10.0mm	84mm
	LCP503208	AG-008L	Locking Plate Radio Distal Dorsal 3.5mm - 8H Left	2.0mm	10.0mm	108mm
	LCP503210	AG-010L	Locking Plate Radio Distal Dorsal 3.5mm - 10H Left	2.0mm	10.0mm	132mm
42			LOCKING PLATE RADIO DISTAL LATERAL 3.5mm			
	LCP502005	AF-005	Locking Plate Radio Distal Lateral 3.5mm - 5H	2.0mm	10.0mm	70mm
	LCP502006	AF-006	Locking Plate Radio Distal Lateral 3.5mm - 6H	2.0mm	10.0mm	84mm
43			LOCKING PLATE RECONSTRUCTION 3.5mm			
	LCP504003	EE-003	Locking Plate Reconstruction 3.5mm - 3H	2.5mm	10.5mm	34mm
	LCP504004	EE-004	Locking Plate Reconstruction 3.5mm - 4H	2.5mm	10.5mm	46mm
	LCP504005	EE-005	Locking Plate Reconstruction 3.5mm - 5H	2.5mm	10.5mm	58mm
	LCP504006	EE-006	Locking Plate Reconstruction 3.5mm - 6H	2.5mm	10.5mm	70mm
	LCP504007	EE-007	Locking Plate Reconstruction 3.5mm - 7H	2.5mm	10.5mm	82mm
	LCP504008	EE-008	Locking Plate Reconstruction 3.5mm - 8H	2.5mm	10.5mm	94mm
	LCP504009	EE-009	Locking Plate Reconstruction 3.5mm - 9H	2.5mm	10.5mm	106mm
	LCP504010	EE-010	Locking Plate Reconstruction 3.5mm - 10H	2.5mm	10.5mm	118mm
	LCP504011	EE-011	Locking Plate Reconstruction 3.5mm - 11H	2.5mm	10.5mm	130mm

	LCP504012	EE-012	Locking Plate Reconstruction 3.5mm - 12H	2.5mm	10.5mm	142mm
	LCP504013	EE-013	Locking Plate Reconstruction 3.5mm - 13H	2.5mm	10.5mm	154mm
	LCP504014	EE-014	Locking Plate Reconstruction 3.5mm - 14H	2.5mm	10.5mm	166mm
	LCP504015	EE-015	Locking Plate Reconstruction 3.5mm - 15H	2.5mm	10.5mm	178mm
	LCP504016	EE-016	Locking Plate Reconstruction 3.5mm - 16H	2.5mm	10.5mm	190mm
44			LOCKING PLATE RECONSTRUCTION 4.5mm			
	LCP555003	AP-003	Locking Plate Reconstruction 4.5mm- 3H	3.0mm	12.0mm	45mm
	LCP555004	AP-004	Locking Plate Reconstruction 4.5mm- 4H	3.0mm	12.0mm	61mm
	LCP555005	AP-005	Locking Plate Reconstruction 4.5mm- 5H	3.0mm	12.0mm	77mm
	LCP555006	AP-006	Locking Plate Reconstruction 4.5mm- 6H	3.0mm	12.0mm	93mm
	LCP555007	AP-007	Locking Plate Reconstruction 4.5mm- 7H	3.0mm	12.0mm	109mm
	LCP555008	AP-008	Locking Plate Reconstruction 4.5mm- 8H	3.0mm	12.0mm	125mm
	LCP555009	AP-009	Locking Plate Reconstruction 4.5mm- 9H	3.0mm	12.0mm	141mm
	LCP555010	AP-010	Locking Plate Reconstruction 4.5mm- 10H	3.0mm	12.0mm	157mm
	LCP555011	AP-011	Locking Plate Reconstruction 4.5mm- 11H	3.0mm	12.0mm	173mm
	LCP555012	AP-012	Locking Plate Reconstruction 4.5mm- 12H	3.0mm	12.0mm	189mm
	LCP555013	AP-013	Locking Plate Reconstruction 4.5mm- 13H	3.0mm	12.0mm	205mm
	LCP555014	AP-014	Locking Plate Reconstruction 4.5mm- 14H	3.0mm	12.0mm	221mm
	LCP555015	AP-015	Locking Plate Reconstruction 4.5mm- 15H	3.0mm	12.0mm	237mm
	LCP555016	AP-016	Locking Plate Reconstruction 4.5mm- 16H	3.0mm	12.0mm	253mm
45			LOCKING PLATE PROXIMAL FEMORAL LATERAL 4.5/5.0mm TYPE-I, RIGHT/LEFT			
	LCP531105	AO-005R	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 5H Right	5.6mm	18.5mm	156mm
	LCP531107	AO-007R	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 7H Right	5.6mm	18.5mm	196mm
	LCP531109	AO-009R	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 9H Right	5.6mm	18.5mm	236mm
	LCP531111	AO-011R	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 11H Right	5.6mm	18.5mm	276mm
	LCP531113	AO-013R	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 13H Right	5.6mm	18.5mm	310mm
	LCP531205	AO-005L	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 5H Left	5.6mm	18.5mm	156mm
	LCP531207	AO-007L	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 7H Left	5.6mm	18.5mm	196mm
	LCP531209	AO-009L	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 9H Left	5.6mm	18.5mm	236mm
	LCP531211	AO-011L	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 11H Left	5.6mm	18.5mm	276mm
	LCP531213	AO-013L	Locking Plate Proximal Femoral Lateral 4.5/5.0mm, Type-I - 13H Left	5.6mm	18.5mm	310mm
46			LOCKING PLATE PROXIMAL FEMORAL HOOK 4.5/5.0mm TYPE-II, RIGHT/LEFT			
	LCP667102	BZ-002R	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 2H Right	6.0mm	17.5mm	133mm
	LCP667104	BZ-004R	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 4H Right	6.0mm	17.5mm	169mm
	LCP667106	BZ-006R	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 6H Right	6.0mm	17.5mm	205mm
	LCP667108	BZ-008R	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 8H Right	6.0mm	17.5mm	241mm
	LCP667110	BZ-010R	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 10H Right	6.0mm	17.5mm	277mm

	LCP667112	BZ-012R	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 12H Right	6.0mm	17.5mm	313mm
	LCP667202	BZ-002L	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 2H Left	6.0mm	17.5mm	133mm
	LCP667204	BZ-004L	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 4H Left	6.0mm	17.5mm	169mm
	LCP667206	BZ-006L	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 6H Left	6.0mm	17.5mm	205mm
	LCP667208	BZ-008L	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 8H Left	6.0mm	17.5mm	241mm
	LCP667210	BZ-010L	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 10H Left	6.0mm	17.5mm	277mm
	LCP667212	BZ-012L	Locking Plate Proximal Femoral Hook 4.5/5.0mm, Type-II - 12H Left	6.0mm	17.5mm	313mm
47			LOCKING PLATE PEDIATRIC HIP 3.5 & 5.0mm			
	LCP6673100	BP-100-75	Locking Plate Pediatric Hip 3.5mm, 100°, Length - 75mm	3.4mm	10.1mm	75mm
	LCP6673110	BP-110-75	Locking Plate Pediatric Hip 3.5mm, 110°, Length - 75mm	3.4mm	10.1mm	75mm
	LCP6673120	BP-120-75	Locking Plate Pediatric Hip 3.5mm, 120°, Length - 75mm	3.4mm	10.1mm	75mm
	LCP6673150	BP-150-60	Locking Plate Pediatric Hip 3.5mm, 150°, Length - 60mm	3.4mm	10.1mm	60mm
	LCP6673140	BP-140-70	Locking Plate Pediatric Hip 3.5mm, 140°, Length - 70mm	3.4mm	10.1mm	70mm
	LCP6674100	BP-100-90	Locking Plate Pediatric Hip 5.0mm, 100°, Length - 90mm	5.5mm	16mm	90mm
	LCP6674110	BP-110-90	Locking Plate Pediatric Hip 5.0mm, 110°, Length - 90mm	5.5mm	16mm	90mm
	LCP6674120	BP-120-95	Locking Plate Pediatric Hip 5.0mm, 120°, Length - 95mm	5.5mm	16mm	95mm
	LCP6674150	BP-150-75	Locking Plate Pediatric Hip 5.0mm, 150°, Length - 75mm	5.5mm	16mm	75mm
	LCP6674140	BP-140-90	Locking Plate Pediatric Hip 5.0mm, 140°, Length - 90mm	5.5mm	16mm	90mm
48			LOCKING PLATE DHS 135° WITH LONG BARREL 38mm			
	LCP410004	AB-004	Locking Plate DHS 135° With Long Barrel - 4H	6.8mm	19.0mm	78mm
	LCP410005	AB-005	Locking Plate DHS 135° With Long Barrel - 5H	6.8mm	19.0mm	94mm
	LCP410006	AB-006	Locking Plate DHS 135° With Long Barrel - 6H	6.8mm	19.0mm	110mm
	LCP410007	AB-007	Locking Plate DHS 135° With Long Barrel - 7H	6.8mm	19.0mm	126mm
	LCP410008	AB-008	Locking Plate DHS 135° With Long Barrel - 8H	6.8mm	19.0mm	142mm
	LCP410009	AB-009	Locking Plate DHS 135° With Long Barrel - 9H	6.8mm	19.0mm	158mm
	LCP410010	AB-010	Locking Plate DHS 135° With Long Barrel - 10H	6.8mm	19.0mm	174mm
	LCP410011	AB-011	Locking Plate DHS 135° With Long Barrel - 11H	6.8mm	19.0mm	190mm
	LCP410012	AB-012	Locking Plate DHS 135° With Long Barrel - 12H	6.8mm	19.0mm	206mm
	LCP410013	AB-013	Locking Plate DHS 135° With Long Barrel - 13H	6.8mm	19.0mm	222mm
	LCP410014	AB-014	Locking Plate DHS 135° With Long Barrel - 14H	6.8mm	19.0mm	238mm
	LCP410015	AB-015	Locking Plate DHS 135° With Long Barrel - 15H	6.8mm	19.0mm	254mm
	LCP410016	AB-016	Locking Plate DHS 135° With Long Barrel - 16H	6.8mm	19.0mm	270mm
	LCP410017	AB-017	Locking Plate DHS 135° With Long Barrel - 17H	6.8mm	19.0mm	286mm
	LCP410018	AB-018	Locking Plate DHS 135° With Long Barrel - 18H	6.8mm	19.0mm	302mm
	LCP410019	AB-019	Locking Plate DHS 135° With Long Barrel - 19H	6.8mm	19.0mm	318mm
	LCP410020	AB-020	Locking Plate DHS 135° With Long Barrel - 20H	6.8mm	19.0mm	334mm
	LCP410021	AB-021	Locking Plate DHS 135° With Long Barrel - 21H	6.8mm	19.0mm	350mm
	LCP410022	AB-022	Locking Plate DHS 135° With Long Barrel - 22H	6.8mm	19.0mm	366mm

49			LOCKING PLATE DCS 95°			
	LCP510006	AC-006	Locking Plate DCS 95° - 6H	6.8mm	16.0mm	100mm
	LCP510007	AC-007	Locking Plate DCS 95° - 7H	6.8mm	16.0mm	115mm
	LCP510008	AC-008	Locking Plate DCS 95° - 8H	6.8mm	16.0mm	130mm
	LCP510009	AC-009	Locking Plate DCS 95° - 9H	6.8mm	16.0mm	146mm
	LCP510010	AC-010	Locking Plate DCS 95° - 10H	6.8mm	16.0mm	163mm
	LCP510011	AC-011	Locking Plate DCS 95° - 11H	6.8mm	16.0mm	180mm
	LCP510012	AC-012	Locking Plate DCS 95° - 12H	6.8mm	16.0mm	198mm
	LCP510013	AC-013	Locking Plate DCS 95° - 13H	6.8mm	16.0mm	211mm
	LCP510014	AC-014	Locking Plate DCS 95° - 14H	6.8mm	16.0mm	225mm
	LCP510015	AC-015	Locking Plate DCS 95° - 15H	6.8mm	16.0mm	242mm
	LCP510016	AC-016	Locking Plate DCS 95° - 16H	6.8mm	16.0mm	260mm
	LCP510017	AC-017	Locking Plate DCS 95° - 17H	6.8mm	16.0mm	283mm
	LCP510018	AC-018	Locking Plate DCS 95° - 18H	6.8mm	16.0mm	306mm
	LCP510019	AC-019	Locking Plate DCS 95° - 19H	6.8mm	16.0mm	322mm
	LCP510020	AC-020	Locking Plate DCS 95° - 20H	6.8mm	16.0mm	338mm
	LCP510021	AC-021	Locking Plate DCS 95° - 21H	6.8mm	16.0mm	354mm
	LCP510022	AC-022	Locking Plate DCS 95° - 22H	6.8mm	16.0mm	370mm
50			LOCKING PLATE CONDYLAR BUTTERESS 4.5/5.0mm, RIGHT/LEFT			
	LCP667506	CA-006R	Locking Plate Condylar Butteress 4.5/5.0mm - 6H Right	6.4mm	11.2mm	170mm
	LCP667508	CA-008R	Locking Plate Condylar Butteress 4.5/5.0mm - 8H Right	6.4mm	11.2mm	206mm
	LCP667510	CA-010R	Locking Plate Condylar Butteress 4.5/5.0mm - 10H Right	6.4mm	11.2mm	242mm
	LCP667512	CA-012R	Locking Plate Condylar Butteress 4.5/5.0mm - 12H Right	6.4mm	11.2mm	278mm
	LCP667514	CA-014R	Locking Plate Condylar Butteress 4.5/5.0mm - 14H Right	6.4mm	11.2mm	314mm
	LCP667606	CA-006L	Locking Plate Condylar Butteress 4.5/5.0mm - 6H Left	6.4mm	11.2mm	170mm
	LCP667608	CA-008L	Locking Plate Condylar Butteress 4.5/5.0mm - 8H Left	6.4mm	11.2mm	206mm
	LCP667610	CA-010L	Locking Plate Condylar Butteress 4.5/5.0mm - 10H Left	6.4mm	11.2mm	242mm
	LCP667612	CA-012L	Locking Plate Condylar Butteress 4.5/5.0mm - 12H Left	6.4mm	11.2mm	278mm
	LCP667614	CA-014L	Locking Plate Condylar Butteress 4.5/5.0mm - 14H Left	6.4mm	11.2mm	314mm
51			LOCKING PLATE NARROW 4.5mm/5.0mm			
	LCP521004	II-004	Locking Plate Narrow 4.5/5.0mm - 4H	4.5mm	13.5mm	70mm
	LCP521005	II-005	Locking Plate Narrow 4.5/5.0mm - 5H	4.5mm	13.5mm	88mm
	LCP521006	II-006	Locking Plate Narrow 4.5/5.0mm - 6H	4.5mm	13.5mm	106mm
	LCP521007	II-007	Locking Plate Narrow 4.5/5.0mm - 7H	4.5mm	13.5mm	124mm
	LCP521008	II-008	Locking Plate Narrow 4.5/5.0mm - 8H	4.5mm	13.5mm	142mm
	LCP521009	II-009	Locking Plate Narrow 4.5/5.0mm - 9H	4.5mm	13.5mm	160mm
	LCP521010	II-010	Locking Plate Narrow 4.5/5.0mm - 10H	4.5mm	13.5mm	178mm
	LCP521011	II-011	Locking Plate Narrow 4.5/5.0mm - 11H	4.5mm	13.5mm	196mm
	LCP521012	II-012	Locking Plate Narrow 4.5/5.0mm - 12H	4.5mm	13.5mm	214mm
	LCP521014	II-014	Locking Plate Narrow 4.5/5.0mm - 14H	4.5mm	13.5mm	250mm

	LCP521016	II-016	Locking Plate Narrow 4.5/5.0mm - 16H	4.5mm	13.5mm	286mm
	LCP521018	II-018	Locking Plate Narrow 4.5/5.0mm - 18H	4.5mm	13.5mm	322mm
	LCP521020	II-020	Locking Plate Narrow 4.5/5.0mm - 20H	4.5mm	13.5mm	358mm
	LCP521022	II-022	Locking Plate Narrow 4.5/5.0mm - 22H	4.5mm	13.5mm	394mm
52			LOCKING PLATE CURVED NARROW 4.5mm/5.0mm			
	LCPC521010	II-010C	Locking Plate Curved Narrow 4.5/5.0mm - 10H	4.5mm	13.5mm	214mm
	LCPC521011	II-011C	Locking Plate Curved Narrow 4.5/5.0mm - 11H	4.5mm	13.5mm	232mm
	LCPC521012	II-012C	Locking Plate Curved Narrow 4.5/5.0mm - 12H	4.5mm	13.5mm	250mm
	LCPC521014	II-014C	Locking Plate Curved Narrow 4.5/5.0mm - 14H	4.5mm	13.5mm	286mm
	LCPC521016	II-016C	Locking Plate Curved Narrow 4.5/5.0mm - 16H	4.5mm	13.5mm	322mm
	LCPC521018	II-018C	Locking Plate Curved Narrow 4.5/5.0mm - 18H	4.5mm	13.5mm	358mm
	LCPC521020	II-020C	Locking Plate Curved Narrow 4.5/5.0mm - 20H	4.5mm	13.5mm	394mm
	LCPC521022	II-022C	Locking Plate Curved Narrow 4.5/5.0mm - 22H	4.5mm	13.5mm	430mm
53			LOCKING PLATE BROAD 4.5mm/5.0mm			
	LCP522005	JJ-005	Locking Plate Broad 4.5/5.0mm - 5H	5.6mm	17.5mm	88mm
	LCP522006	JJ-006	Locking Plate Broad 4.5/5.0mm - 6H	5.6mm	17.5mm	106mm
	LCP522007	JJ-007	Locking Plate Broad 4.5/5.0mm - 7H	5.6mm	17.5mm	124mm
	LCP522008	JJ-008	Locking Plate Broad 4.5/5.0mm - 8H	5.6mm	17.5mm	142mm
	LCP522009	JJ-009	Locking Plate Broad 4.5/5.0mm - 9H	5.6mm	17.5mm	160mm
	LCP522010	JJ-010	Locking Plate Broad 4.5/5.0mm - 10H	5.6mm	17.5mm	178mm
	LCP522011	JJ-011	Locking Plate Broad 4.5/5.0mm - 11H	5.6mm	17.5mm	196mm
	LCP522012	JJ-012	Locking Plate Broad 4.5/5.0mm - 12H	5.6mm	17.5mm	214mm
	LCP522014	JJ-014	Locking Plate Broad 4.5/5.0mm - 14H	5.6mm	17.5mm	250mm
	LCP522016	JJ-016	Locking Plate Broad 4.5/5.0mm - 16H	5.6mm	17.5mm	286mm
	LCP522018	JJ-018	Locking Plate Broad 4.5/5.0mm - 18H	5.6mm	17.5mm	322mm
	LCP522020	JJ-020	Locking Plate Broad 4.5/5.0mm - 20H	5.6mm	17.5mm	358mm
	LCP522022	JJ-022	Locking Plate Broad 4.5/5.0mm - 22H	5.6mm	17.5mm	394mm
54			LOCKING PLATE CURVED BROAD 4.5mm/5.0mm			
	LCPC522010	JJ-010C	Locking Plate Curved Broad 4.5/5.0mm - 10H	5.6mm	17.5mm	214mm
	LCPC522011	JJ-011C	Locking Plate Curved Broad 4.5/5.0mm - 11H	5.6mm	17.5mm	250mm
	LCPC522012	JJ-012C	Locking Plate Curved Broad 4.5/5.0mm - 12H	5.6mm	17.5mm	286mm
	LCPC522014	JJ-014C	Locking Plate Curved Broad 4.5/5.0mm - 14H	5.6mm	17.5mm	322mm
	LCPC522016	JJ-016C	Locking Plate Curved Broad 4.5/5.0mm - 16H	5.6mm	17.5mm	358mm
	LCPC522018	JJ-018C	Locking Plate Curved Broad 4.5/5.0mm - 18H	5.6mm	17.5mm	394mm
	LCPC522020	JJ-020C	Locking Plate Curved Broad 4.5/5.0mm - 20H	5.6mm	17.5mm	430mm
	LCPC522022	JJ-022C	Locking Plate Curved Broad 4.5/5.0mm - 22H	5.6mm	17.5mm	466mm
	LCPC522024	JJ-024C	Locking Plate Curved Broad 4.5/5.0mm - 24H	5.6mm	17.5mm	502mm
	LCPC522026	JJ-026C	Locking Plate Curved Broad 4.5/5.0mm - 26H	5.6mm	17.5mm	538mm
55			LOCKING PLATE DISTAL FEMORAL 4.5/5.0mm, RIGHT/LEFT			

	LCP527105	KK-005	Locking Plate Distal Femoral 4.5/5.0mm - 5H Right	6.5mm	16.5mm	156mm
	LCP527107	KK-007	Locking Plate Distal Femoral 4.5/5.0mm - 7H Right	6.5mm	16.5mm	196mm
	LCP527109	KK-009	Locking Plate Distal Femoral 4.5/5.0mm - 9H Right	6.5mm	16.5mm	236mm
	LCP527111	KK-011	Locking Plate Distal Femoral 4.5/5.0mm - 11H Right	6.5mm	16.5mm	276mm
	LCP527113	KK-013	Locking Plate Distal Femoral 4.5/5.0mm - 13H Right	6.5mm	16.5mm	316mm
	LCP527205	LL-005	Locking Plate Distal Femoral 4.5/5.0mm - 5H Left	6.5mm	16.5mm	156mm
	LCP527207	LL-007	Locking Plate Distal Femoral 4.5/5.0mm - 7H Left	6.5mm	16.5mm	196mm
	LCP527209	LL-009	Locking Plate Distal Femoral 4.5/5.0mm - 9H Left	6.5mm	16.5mm	236mm
	LCP527211	LL-011	Locking Plate Distal Femoral 4.5/5.0mm - 11H Left	6.5mm	16.5mm	276mm
	LCP527213	LL-013	Locking Plate Distal Femoral 4.5/5.0mm - 13H Left	6.5mm	16.5mm	316mm
56			LOCKING PLATE T-FIX DISTAL FEMORAL 4.5/5.0mm, RIGHT/LEFT			
	LCP531104	TFD-004R	Locking Plate T-Fix Distal Femoral 4.5/5.0mm - 4H, Right	6.5mm	16.5mm	121mm
	LCP531204	TFD-004L	Locking Plate T-Fix Distal Femoral 4.5/5.0mm - 4H, Left	6.5mm	16.5mm	121mm
57			Locking Plate Periarticular Distal Femur 4.5/5.0mm Left & Right			
	LCP527306	PADF-006R	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 6H Right	6.5mm	16.5mm	196mm
	LCP527308	PADF-008R	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 8H Right	6.5mm	16.5mm	236mm
	LCP527310	PADF-010R	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 10H Right	6.5mm	16.5mm	276mm
	LCP527312	PADF-012R	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 12H Right	6.5mm	16.5mm	316mm
	LCP527314	PADF-014R	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 14H Right	6.5mm	16.5mm	356mm
	LCP527316	PADF-016R	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 16H Right	6.5mm	16.5mm	396mm
	LCP527406	PADF-006L	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 6H Left	6.5mm	16.5mm	196mm
	LCP527408	PADF-008L	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 8H Left	6.5mm	16.5mm	236mm
	LCP527410	PADF-010L	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 10H Left	6.5mm	16.5mm	276mm
	LCP527412	PADF-012L	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 11H Left	6.5mm	16.5mm	316mm
	LCP527414	PADF-014L	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 13H Left	6.5mm	16.5mm	356mm
	LCP527416	PADF-016L	Locking Plate Periarticular Distal Femur 4.5/5.0mm - 15H Left	6.5mm	16.5mm	396mm
58			LOCKING T- PLATE SINGLE CURVED 4.5/5.0mm			
	LCP524003	UU-003	Locking T - Plate Single Curved - 3H	2.5mm	12.4mm	68mm
	LCP524004	UU-004	Locking T - Plate Single Curved - 4H	2.5mm	12.4mm	84mm
	LCP524005	UU-005	Locking T - Plate Single Curved - 5H	2.5mm	12.4mm	100mm
	LCP524006	UU-006	Locking T - Plate Single Curved - 6H	2.5mm	12.4mm	116mm
	LCP524007	UU-007	Locking T - Plate Single Curved - 7H	2.5mm	12.4mm	132mm
	LCP524008	UU-008	Locking T - Plate Single Curved - 8H	2.5mm	12.4mm	148mm
	LCP524010	UU-010	Locking T - Plate Single Curved - 10H	2.5mm	12.4mm	180mm
59			LOCKING T - PLATE DOUBLE CURVED 4.5/5.0mm			
	LCP523003	VV-003	Locking T-Plate Double Curved - 3H	3.3mm	12.4mm	68mm
	LCP523004	VV-004	Locking T-Plate Double Curved - 4H	3.3mm	12.4mm	84mm

	LCP523005	VV-005	Locking T-Plate Double Curved - 5H	3.3mm	12.4mm	100mm
	LCP523006	VV-006	Locking T-Plate Double Curved - 6H	3.3mm	12.4mm	116mm
	LCP523007	VV-007	Locking T-Plate Double Curved - 7H	3.3mm	12.4mm	132mm
	LCP523008	VV-008	Locking T-Plate Double Curved - 8H	3.3mm	12.4mm	148mm
	LCP523010	VV-010	Locking T-Plate Double Curved - 10H	3.3mm	12.4mm	164mm
60			LOCKING T-FIX T-PLATE 4.5/5.0mm			
	LCP537004	TFT-004	Locking T-Fix T-Plate 4.5/5.0mm - 4H	2.5mm	17.0mm	115mm
61			LOCKING L -PLATE DOUBLE CURVED 4.5/ 5.0mm, RIGHT/LEFT			
	LCP525103	WX-003	Locking L- Plate Double Curved 4.5/5.0mm - 3 H Right	2.5mm	15.6mm	69mm
	LCP525104	WX-004	Locking L- Plate Double Curved 4.5/5.0mm - 4 H Right	2.5mm	15.6mm	85mm
	LCP525105	WX-005	Locking L- Plate Double Curved 4.5/5.0mm - 5 H Right	2.5mm	15.6mm	101mm
	LCP525106	WX-006	Locking L- Plate Double Curved 4.5/5.0mm - 6 H Right	2.5mm	15.6mm	117mm
	LCP525107	WX-007	Locking L- Plate Double Curved 4.5/5.0mm - 7 H Right	2.5mm	15.6mm	133mm
	LCP525108	WX-008	Locking L- Plate Double Curved 4.5/5.0mm - 8 H Right	2.5mm	15.6mm	149mm
	LCP525203	WW-003	Locking L- Plate Double Curved 4.5/5.0mm - 3 H Left	2.5mm	15.6mm	69mm
	LCP525204	WW-004	Locking L- Plate Double Curved 4.5/5.0mm - 4 H Left	2.5mm	15.6mm	85mm
	LCP525205	WW-005	Locking L- Plate Double Curved 4.5/5.0mm - 5 H Left	2.5mm	15.6mm	101mm
	LCP525206	WW-006	Locking L- Plate Double Curved 4.5/5.0mm - 6 H Left	2.5mm	15.6mm	117mm
	LCP525207	WW-007	Locking L- Plate Double Curved 4.5/5.0mm - 7 H Left	2.5mm	15.6mm	133mm
	LCP525208	WW-008	Locking L- Plate Double Curved 4.5/5.0mm - 8 H Left	2.5mm	15.6mm	149mm
62			LOCKING PLATE PROXIMAL LATERAL TIBIA 4.5/5.0mm, RIGHT/LEFT			
	LCP526105	MM-005	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 5H, Right	4.5mm	11.2mm	140mm
	LCP526107	MM-007	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 7H, Right	4.5mm	11.2mm	180mm
	LCP526109	MM-009	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 9H, Right	4.5mm	11.2mm	220mm
	LCP526111	MM-011	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 11H, Right	4.5mm	11.2mm	260mm
	LCP526113	MM-013	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 13H, Right	4.5mm	11.2mm	300mm
	LCP526205	NN-005	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 5H, Left	4.5mm	11.2mm	140mm
	LCP526207	NN-007	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 7H, Left	4.5mm	11.2mm	180mm
	LCP526209	NN-009	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 9H, Left	4.5mm	11.2mm	220mm
	LCP526211	NN-011	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 11H, Left	4.5mm	11.2mm	260mm
	LCP526213	NN-013	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 13H, Left	4.5mm	11.2mm	300mm
63			LOCKING PLATE T-FIX PROXIMAL TIBIA 4.5/5.0mm, RIGHT/LEFT			
	LCP530103	TFP-003R	Locking Plate T-Fix Proximal Tibia 4.5/5.0mm - 3H, Right	4.5mm	11.2mm	102mm
	LCP530203	TFP-003L	Locking Plate T-Fix Proximal Tibia 4.5/5.0mm - 3H, Left	4.5mm	11.2mm	102mm
64			LOCKING PLATE PROXIMAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP667704	BQ-004R	Locking Plate Proximal Tibia 3.5mm - 4 H Right	4.0mm	10.7mm	81mm

	LCP667706	BQ-006R	Locking Plate Proximal Tibia 3.5mm - 6 H Right	4.0mm	10.7mm	107mm
	LCP667708	BQ-008R	Locking Plate Proximal Tibia 3.5mm - 8 H Right	4.0mm	10.7mm	133mm
	LCP667710	BQ-010R	Locking Plate Proximal Tibia 3.5mm - 10 H Right	4.0mm	10.7mm	159mm
	LCP667712	BQ-012R	Locking Plate Proximal Tibia 3.5mm - 12 H Right	4.0mm	10.7mm	185mm
	LCP667714	BQ-014R	Locking Plate Proximal Tibia 3.5mm - 14 H Right	4.0mm	10.7mm	211mm
	LCP667716	BQ-016R	Locking Plate Proximal Tibia 3.5mm - 16 H Right	4.0mm	10.7mm	237mm
	LCP667804	BQ-004L	Locking Plate Proximal Tibia 3.5mm - 4 H Left	4.0mm	10.7mm	81mm
	LCP667806	BQ-006L	Locking Plate Proximal Tibia 3.5mm - 6 H Left	4.0mm	10.7mm	107mm
	LCP667808	BQ-008L	Locking Plate Proximal Tibia 3.5mm - 8 H Left	4.0mm	10.7mm	133mm
	LCP667810	BQ-010L	Locking Plate Proximal Tibia 3.5mm - 10 H Left	4.0mm	10.7mm	159mm
	LCP667812	BQ-012L	Locking Plate Proximal Tibia 3.5mm - 12 H Left	4.0mm	10.7mm	185mm
	LCP667814	BQ-014L	Locking Plate Proximal Tibia 3.5mm - 14 H Left	4.0mm	10.7mm	211mm
	LCP667816	BQ-016L	Locking Plate Proximal Tibia 3.5mm - 16 H Left	4.0mm	10.7mm	237mm
65			LOCKING PLATE PROXIMAL LATERAL TIBIA 4.5/5.0mm, RIGHT/LEFT			
	LCP667904	CB-004R	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 4 H Right	4.0mm	13.0mm	82mm
	LCP667906	CB-006R	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 6 H Right	4.0mm	13.0mm	118mm
	LCP667908	CB-008R	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 8 H Right	4.0mm	13.0mm	154mm
	LCP667910	CB-010R	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 10 H Right	4.0mm	13.0mm	190mm
	LCP667912	CB-012R	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 12 H Right	4.0mm	13.0mm	226mm
	LCP668004	CB-004L	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 4 H Left	4.0mm	13.0mm	82mm
	LCP668006	CB-006L	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 6 H Left	4.0mm	13.0mm	118mm
	LCP668008	CB-008L	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 8 H Left	4.0mm	13.0mm	154mm
	LCP668010	CB-010L	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 10 H Left	4.0mm	13.0mm	190mm
	LCP668012	CB-012L	Locking Plate Proximal Lateral Tibia 4.5/5.0mm - 12 H Left	4.0mm	13.0mm	226mm
66			LOCKING PLATE MEDIAL PROXIMAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP881004	BR-004R	Locking Plate Medial Proximal Tibia 3.5mm - 4 H Right	4.0mm	11.0mm	93mm
	LCP881006	BR-006R	Locking Plate Medial Proximal Tibia 3.5mm - 6 H Right	4.0mm	11.0mm	119mm
	LCP881008	BR-008R	Locking Plate Medial Proximal Tibia 3.5mm - 8 H Right	4.0mm	11.0mm	145mm
	LCP881010	BR-010R	Locking Plate Medial Proximal Tibia 3.5mm - 10 H Right	4.0mm	11.0mm	171mm
	LCP881012	BR-012R	Locking Plate Medial Proximal Tibia 3.5mm - 12 H Right	4.0mm	11.0mm	197mm
	LCP881014	BR-014R	Locking Plate Medial Proximal Tibia 3.5mm - 14 H Right	4.0mm	11.0mm	223mm
	LCP880104	BR-004L	Locking Plate Medial Proximal Tibia 3.5mm - 4 H Left	4.0mm	11.0mm	93mm
	LCP880106	BR-006L	Locking Plate Medial Proximal Tibia 3.5mm - 6 H Left	4.0mm	11.0mm	119mm
	LCP880108	BR-008L	Locking Plate Medial Proximal Tibia 3.5mm - 8 H Left	4.0mm	11.0mm	145mm
	LCP880110	BR-010L	Locking Plate Medial Proximal Tibia 3.5mm - 10 H Left	4.0mm	11.0mm	171mm
	LCP880112	BR-012L	Locking Plate Medial Proximal Tibia 3.5mm - 12 H Left	4.0mm	11.0mm	197mm
	LCP880114	BR-014L	Locking Plate Medial Proximal Tibia 3.5mm - 14 H Left	4.0mm	11.0mm	223mm

67			LOCKING PLATE MEDIAL PROXIMAL TIBIA 3.5mm TYPE-II LOW BEND, RIGHT/LEF			
	LCP681304	DD-004R	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 4H Right	3.8mm	10.7mm	81mm
	LCP681306	DD-006R	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 6H Right	3.8mm	10.7mm	107mm
	LCP681308	DD-008R	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 8H Right	3.8mm	10.7mm	133mm
	LCP681310	DD-010R	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 10H Right	3.8mm	10.7mm	159mm
	LCP681312	DD-012R	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 12H Right	3.8mm	10.7mm	185mm
	LCP681314	DD-014R	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 14H Right	3.8mm	10.7mm	211mm
	LCP681404	DD-004L	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 4H Left	3.8mm	10.7mm	81mm
	LCP681406	DD-006L	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 6H Left	3.8mm	10.7mm	107mm
	LCP681408	DD-008L	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 8H Left	3.8mm	10.7mm	133mm
	LCP681410	DD-010L	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 10H Left	3.8mm	10.7mm	159mm
	LCP681412	DD-012L	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 12H Left	3.8mm	10.7mm	185mm
	LCP681414	DD-014L	Locking Plate Proximal Tibia 3.5mm, Type-II Low Bend - 14H Left	3.8mm	10.7mm	211mm
68			LOCKING PLATE PERIARTICULAR PROXIMAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP691104	PAPT-3504R	Locking Plate Periarticular Proximal Tibia 3.5mm - 4 H Right	3.8mm	11.0mm	80mm
	LCP691105	PAPT-3505R	Locking Plate Periarticular Proximal Tibia 3.5mm - 5 H Right	3.8mm	11.0mm	92mm
	LCP691106	PAPT-3506R	Locking Plate Periarticular Proximal Tibia 3.5mm - 6 H Right	3.8mm	11.0mm	104mm
	LCP691107	PAPT-3507R	Locking Plate Periarticular Proximal Tibia 3.5mm - 7 H Right	3.8mm	11.0mm	116mm
	LCP691108	PAPT-3508R	Locking Plate Periarticular Proximal Tibia 3.5mm - 8 H Right	3.8mm	11.0mm	128mm
	LCP691109	PAPT-3509R	Locking Plate Periarticular Proximal Tibia 3.5mm - 9 H Right	3.8mm	11.0mm	140mm
	LCP691110	PAPT-3510R	Locking Plate Periarticular Proximal Tibia 3.5mm - 10 H Right	3.8mm	11.0mm	152mm
	LCP691111	PAPT-3511R	Locking Plate Periarticular Proximal Tibia 3.5mm - 11 H Right	3.8mm	11.0mm	164mm
	LCP691112	PAPT-3512R	Locking Plate Periarticular Proximal Tibia 3.5mm - 12 H Right	3.8mm	11.0mm	176mm
	LCP691113	PAPT-3513R	Locking Plate Periarticular Proximal Tibia 3.5mm - 13 H Right	3.8mm	11.0mm	188mm
	LCP691114	PAPT-3514R	Locking Plate Periarticular Proximal Tibia 3.5mm - 14 H Right	3.8mm	11.0mm	200mm
	LCP691204	PAPT-3504L	Locking Plate Periarticular Proximal Tibia 3.5mm - 4 H Left	3.8mm	11.0mm	80mm
	LCP691205	PAPT-3505L	Locking Plate Periarticular Proximal Tibia 3.5mm - 5 H Left	3.8mm	11.0mm	92mm
	LCP691206	PAPT-3506L	Locking Plate Periarticular Proximal Tibia 3.5mm - 6 H Left	3.8mm	11.0mm	104mm
	LCP691207	PAPT-3507L	Locking Plate Periarticular Proximal Tibia 3.5mm - 7 H Left	3.8mm	11.0mm	116mm
	LCP691208	PAPT-3508L	Locking Plate Periarticular Proximal Tibia 3.5mm - 8 H Left	3.8mm	11.0mm	128mm
	LCP691209	PAPT-3509L	Locking Plate Periarticular Proximal Tibia 3.5mm - 9 H Left	3.8mm	11.0mm	140mm
	LCP691210	PAPT-3510L	Locking Plate Periarticular Proximal Tibia 3.5mm - 10 H Left	3.8mm	11.0mm	152mm
	LCP691211	PAPT-3511L	Locking Plate Periarticular Proximal Tibia 3.5mm - 11 H Left	3.8mm	11.0mm	164mm
	LCP691212	PAPT-3512L	Locking Plate Periarticular Proximal Tibia 3.5mm - 12 H Left	3.8mm	11.0mm	176mm
	LCP691213	PAPT-3513L	Locking Plate Periarticular Proximal Tibia 3.5mm - 13 H Left	3.8mm	11.0mm	188mm
	LCP691214	PAPT-3514L	Locking Plate Periarticular Proximal Tibia 3.5mm - 14 H Left	3.8mm	11.0mm	200mm
69			LOCKING PLATE PERIARTICULAR PROXIMAL TIBIA 3.5/5.0mm, RIGHT/LEFT			
	LCP691304	PAPT-5004R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 4 H Right	4.0mm	13.0mm	104mm
	LCP691305	PAPT-5005R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 5 H Right	4.0mm	13.0mm	116mm

	LCP691306	PAPT-5006R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 6 H Right	4.0mm	13.0mm	128mm
	LCP691307	PAPT-5007R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 7 H Right	4.0mm	13.0mm	140mm
	LCP691308	PAPT-5008R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 8 H Right	4.0mm	13.0mm	152mm
	LCP691309	PAPT-5009R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 9 H Right	4.0mm	13.0mm	164mm
	LCP691310	PAPT-5010R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 10 H Right	4.0mm	13.0mm	176mm
	LCP691311	PAPT-5011R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 11 H Right	4.0mm	13.0mm	188mm
	LCP691312	PAPT-5012R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 12 H Right	4.0mm	13.0mm	200mm
	LCP691313	PAPT-5013R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 13 H Right	4.0mm	13.0mm	212mm
	LCP691314	PAPT-5014R	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 14 H Right	4.0mm	13.0mm	220mm
	LCP691404	PAPT-5004L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 4 H Left	4.0mm	13.0mm	104mm
	LCP691405	PAPT-5005L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 5 H Left	4.0mm	13.0mm	116mm
	LCP691406	PAPT-5006L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 6 H Left	4.0mm	13.0mm	128mm
	LCP691407	PAPT-5007L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 7 H Left	4.0mm	13.0mm	140mm
	LCP691408	PAPT-5008L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 8 H Left	4.0mm	13.0mm	152mm
	LCP691409	PAPT-5009L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 9 H Left	4.0mm	13.0mm	164mm
	LCP691410	PAPT-5010L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 10 H Left	4.0mm	13.0mm	176mm
	LCP691411	PAPT-5011L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 11 H Left	4.0mm	13.0mm	188mm
	LCP691412	PAPT-5012L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 12 H Left	4.0mm	13.0mm	200mm
	LCP691413	PAPT-5013L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 13 H Left	4.0mm	13.0mm	212mm
	LCP691414	PAPT-5014L	Locking Plate Periarticular Proximal Tibia 3.5/5.0mm - 14 H Left	4.0mm	13.0mm	220mm
70			LOCKING POSTEROMEDIAL PROXIMAL TIBIA PLATE 3.5mm			
	LCP680001	PPTP-001	Locking Posteromedial Proximal Tibia Plate - 1 Holes	3.8mm	11.0mm	69mm
	LCP680002	PPTP-002	Locking Posteromedial Proximal Tibia Plate - 2 Holes	3.8mm	11.0mm	79mm
	LCP680004	PPTP-004	Locking Posteromedial Proximal Tibia Plate - 4 Holes	3.8mm	11.0mm	105mm
	LCP680006	PPTP-006	Locking Posteromedial Proximal Tibia Plate - 6 Holes	3.8mm	11.0mm	131mm
	LCP680008	PPTP-008	Locking Posteromedial Proximal Tibia Plate - 8 Holes	3.8mm	11.0mm	157mm
	LCP680010	PPTP-010	Locking Posteromedial Proximal Tibia Plate - 10 Holes	3.8mm	11.0mm	183mm
71			LOCKING PLATE MEDIAL DISTAL TIBIA WITHOUT TAB 3.5mm, RIGHT/LEFT			
	LCP668604	BX-004R	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 4H Right	4.0mm	13.2mm	116mm
	LCP668606	BX-006R	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 6H Right	4.0mm	13.2mm	142mm
	LCP668608	BX-008R	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 8H Right	4.0mm	13.2mm	168mm
	LCP668610	BX-010R	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 10H Right	4.0mm	13.2mm	194mm
	LCP668612	BX-012R	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 12H Right	4.0mm	13.2mm	220mm
	LCP668704	BX-004L	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 4H Left	4.0mm	13.2mm	116mm
	LCP668706	BX-006L	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 6H Left	4.0mm	13.2mm	142mm
	LCP668708	BX-008L	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 8H Left	4.0mm	13.2mm	168mm
	LCP668710	BX-010L	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 10H Left	4.0mm	13.2mm	194mm
	LCP668712	BX-012L	Locking Plate Medial Distal Tibia Without Tab 3.5mm - 12H Left	4.0mm	13.2mm	220mm

72			LOCKING PLATE MEDIAL PROXIMAL TIBIA 4.5mm, RIGHT/LEFT			
	LCP681704	MPT-004R	Locking Plate Medial Proximal Tibia 4.5mm - 4 H Right	4.5mm	14.5mm	142mm
	LCP681706	MPT-006R	Locking Plate Medial Proximal Tibia 4.5mm - 6 H Right	4.5mm	14.5mm	168mm
	LCP681708	MPT-008R	Locking Plate Medial Proximal Tibia 4.5mm - 8 H Right	4.5mm	14.5mm	194mm
	LCP681710	MPT-010R	Locking Plate Medial Proximal Tibia 4.5mm - 10 H Right	4.5mm	14.5mm	220mm
	LCP681712	MPT-012R	Locking Plate Medial Proximal Tibia 4.5mm - 12 H Right	4.5mm	14.5mm	246mm
	LCP681804	MPT-004L	Locking Plate Medial Proximal Tibia 4.5mm - 4 H Left	4.5mm	14.5mm	142mm
	LCP681806	MPT-006L	Locking Plate Medial Proximal Tibia 4.5mm - 6 H Left	4.5mm	14.5mm	168mm
	LCP681808	MPT-008L	Locking Plate Medial Proximal Tibia 4.5mm - 8 H Left	4.5mm	14.5mm	194mm
	LCP681810	MPT-010L	Locking Plate Medial Proximal Tibia 4.5mm - 10 H Left	4.5mm	14.5mm	220mm
	LCP681812	MPT-012L	Locking Plate Medial Proximal Tibia 4.5mm - 12 H Left	4.5mm	14.5mm	246mm
73			LOCKING PLATE ANATOMICAL PROXIMAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP677705	APT-005R	Locking Plate Anatomical Proximal Tibia 3.5mm - 5 H Right	4.0mm	13.2mm	116mm
	LCP677707	APT-007R	Locking Plate Anatomical Proximal Tibia 3.5mm - 7 H Right	4.0mm	13.2mm	142mm
	LCP677709	APT-009R	Locking Plate Anatomical Proximal Tibia 3.5mm - 9 H Right	4.0mm	13.2mm	168mm
	LCP677711	APT-011R	Locking Plate Anatomical Proximal Tibia 3.5mm - 11 H Right	4.0mm	13.2mm	194mm
	LCP677713	APT-013R	Locking Plate Anatomical Proximal Tibia 3.5mm - 13 H Right	4.0mm	13.2mm	220mm
	LCP677805	APT-005L	Locking Plate Anatomical Proximal Tibia 3.5mm - 5 H Right	4.0mm	13.2mm	116mm
	LCP677807	APT-007L	Locking Plate Anatomical Proximal Tibia 3.5mm - 7 H Left	4.0mm	13.2mm	142mm
	LCP677809	APT-009L	Locking Plate Anatomical Proximal Tibia 3.5mm - 9 H Left	4.0mm	13.2mm	168mm
	LCP677811	APT-011L	Locking Plate Anatomical Proximal Tibia 3.5mm - 11 H Left	4.0mm	13.2mm	194mm
	LCP677813	APT-013L	Locking Plate Anatomical Proximal Tibia 3.5mm - 13 H Left	4.0mm	13.2mm	220mm
74			LOCKING PLATE ANATOMICAL PROXIMAL TIBIA 4.5/5.0mm, RIGHT/LEFT			
	LCP687105	AR-005R	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 5 H Right	4.5mm	14.5mm	142mm
	LCP687107	AR-007R	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 7 H Right	4.5mm	14.5mm	168mm
	LCP687109	AR-009R	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 9 H Right	4.5mm	14.5mm	194mm
	LCP687111	AR-011R	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 11 H Right	4.5mm	14.5mm	220mm
	LCP687113	AR-013R	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 13 H Right	4.5mm	14.5mm	246mm
	LCP687205	AR-005L	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 5 H Left	4.5mm	14.5mm	142mm
	LCP687207	AR-007L	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 7 H Left	4.5mm	14.5mm	168mm
	LCP687209	AR-009L	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 9 H Left	4.5mm	14.5mm	194mm
	LCP687211	AR-011L	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 11 H Left	4.5mm	14.5mm	220mm
	LCP687213	AR-013L	Locking Plate Anatomical Proximal Tibia 4.5/5.0mm - 13 H Left	4.5mm	14.5mm	246mm
75			LOCKING PLATE DISTAL TIBIA 2.7/3.5mm, RIGHT/LEFT			
	LCP528404	BG-004R	Locking Plate Distal Tibia 2.7/3.5mm - 4H, Right	3.7mm	11.4mm	117mm
	LCP528406	BG-006R	Locking Plate Distal Tibia 2.7/3.5mm - 6H, Right	3.7mm	11.4mm	144mm
	LCP528408	BG-008R	Locking Plate Distal Tibia 2.7/3.5mm - 8H, Right	3.7mm	11.4mm	171mm

	LCP528410	BG-010R	Locking Plate Distal Tibia2.7/3.5mm - 10H, Right	3.7mm	11.4mm	198mm
	LCP528412	BG-012R	Locking Plate Distal Tibia2.7/3.5mm - 12H, Right	3.7mm	11.4mm	225mm
	LCP528304	BG-004L	Locking Plate Distal Tibia2.7/3.5mm - 4H, Left	3.7mm	11.4mm	117mm
	LCP528306	BG-006L	Locking Plate Distal Tibia2.7/3.5mm - 6H, Left	3.7mm	11.4mm	144mm
	LCP528308	BG-008L	Locking Plate Distal Tibia2.7/3.5mm - 8H, Left	3.7mm	11.4mm	171mm
	LCP528310	BG-010L	Locking Plate Distal Tibia2.7/3.5mm - 10H, Left	3.7mm	11.4mm	198mm
	LCP528312	BG-012L	Locking Plate Distal Tibia2.7/3.5mm - 12H, Left	3.7mm	11.4mm	225mm
76			LOCKING PLATE DISTAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP528104	AN-004R	Locking Plate Distal Tibia 3.5mm - 4H, Right	3.7mm	11.4mm	117mm
	LCP528106	AN-006R	Locking Plate Distal Tibia 3.5mm - 6H, Right	3.7mm	11.4mm	144mm
	LCP528108	AN-008R	Locking Plate Distal Tibia 3.5mm - 8H, Right	3.7mm	11.4mm	171mm
	LCP528110	AN-010R	Locking Plate Distal Tibia 3.5mm - 10H, Right	3.7mm	11.4mm	198mm
	LCP528112	AN-012R	Locking Plate Distal Tibia 3.5mm - 12H, Right	3.7mm	11.4mm	225mm
	LCP528204	AN-004L	Locking Plate Distal Tibia 3.5mm - 4H, Left	3.7mm	11.4mm	117mm
	LCP528206	AN-006L	Locking Plate Distal Tibia 3.5mm - 6H, Left	3.7mm	11.4mm	144mm
	LCP528208	AN-008L	Locking Plate Distal Tibia 3.5mm - 8H, Left	3.7mm	11.4mm	171mm
	LCP528210	AN-010L	Locking Plate Distal Tibia 3.5mm - 10H, Left	3.7mm	11.4mm	198mm
	LCP528212	AN-012L	Locking Plate Distal Tibia 3.5mm - 12H, Left	3.7mm	11.4mm	225mm
77			LOCKING PLATE METAPHYSEAL 3.5/4.5/5.0mm			
	LCP507003	AE-003	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 3H	4.0mm	13.5mm	118mm
	LCP507004	AE-004	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 4H	4.0mm	13.5mm	136mm
	LCP507005	AE-005	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 5H	4.0mm	13.5mm	154mm
	LCP507006	AE-006	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 6H	4.0mm	13.5mm	172mm
	LCP507007	AE-007	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 7H	4.0mm	13.5mm	190mm
	LCP507008	AE-008	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 8H	4.0mm	13.5mm	208mm
	LCP507009	AE-009	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 9H	4.0mm	13.5mm	226mm
	LCP507011	AE-011	Locking Plate Metaphyseal 3.5/4.5/5.0mm - 11H	4.0mm	13.5mm	262mm
78			LOCKING PLATE METAPHYSEAL FOR DISTAL MEDIAL TIBIA 3.5/4.5/5.0mm, RIGHT/LEFT			
	LCP678504	BY-004R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 4H Right	4.5mm	15.5mm	124mm
	LCP678505	BY-005R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 5H Right	4.5mm	15.5mm	142mm
	LCP678506	BY-006R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 6H Right	4.5mm	15.5mm	160mm
	LCP678507	BY-007R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 7H Right	4.5mm	15.5mm	178mm
	LCP678508	BY-008R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 8H Right	4.5mm	15.5mm	196mm
	LCP678509	BY-009R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 9H Right	4.5mm	15.5mm	214mm
	LCP678510	BY-010R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 10H Right	4.5mm	15.5mm	232mm
	LCP678512	BY-012R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 12H Right	4.5mm	15.5mm	268mm
	LCP678514	BY-014R	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 14H Right	4.5mm	15.5mm	304mm

	LCP678604	BY-004L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 4H Left	4.5mm	15.5mm	124mm
	LCP678605	BY-005L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 5H Left	4.5mm	15.5mm	142mm
	LCP678606	BY-006L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 6H Left	4.5mm	15.5mm	160mm
	LCP678607	BY-007L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 7H Left	4.5mm	15.5mm	178mm
	LCP678608	BY-008L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 8H Left	4.5mm	15.5mm	196mm
	LCP678609	BY-009L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 9H Left	4.5mm	15.5mm	214mm
	LCP678610	BY-010L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 10H Left	4.5mm	15.5mm	232mm
	LCP678612	BY-012L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 12H Left	4.5mm	15.5mm	268mm
	LCP678614	BY-014L	Locking Plate Metaphyseal For Distal Medial Tibia 3.5/4.5/5.0mm - 14H Left	4.5mm	15.5mm	304mm
79			LOCKING PLATE ANTEROLATERAL DISTAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP668405	BW-005R	Locking Plate Anterolateral Distal Tibia 3.5mm - 5 H Right	4.0mm	13.4mm	80mm
	LCP668407	BW-007R	Locking Plate Anterolateral Distal Tibia 3.5mm - 7 H Right	4.0mm	13.4mm	106mm
	LCP668409	BW-009R	Locking Plate Anterolateral Distal Tibia 3.5mm - 9 H Right	4.0mm	13.4mm	132mm
	LCP668411	BW-011R	Locking Plate Anterolateral Distal Tibia 3.5mm - 11 H Right	4.0mm	13.4mm	158mm
	LCP668413	BW-013R	Locking Plate Anterolateral Distal Tibia 3.5mm - 13 H Right	4.0mm	13.4mm	184mm
	LCP668415	BW-015R	Locking Plate Anterolateral Distal Tibia 3.5mm - 15 H Right	4.0mm	13.4mm	210mm
	LCP668505	BW-005L	Locking Plate Anterolateral Distal Tibia 3.5mm - 5 H Left	4.0mm	13.4mm	80mm
	LCP668507	BW-007L	Locking Plate Anterolateral Distal Tibia 3.5mm - 7 H Left	4.0mm	13.4mm	106mm
	LCP668509	BW-009L	Locking Plate Anterolateral Distal Tibia 3.5mm - 9 H Left	4.0mm	13.4mm	132mm
	LCP668511	BW-011L	Locking Plate Anterolateral Distal Tibia 3.5mm - 11 H Left	4.0mm	13.4mm	158mm
	LCP668513	BW-013L	Locking Plate Anterolateral Distal Tibia 3.5mm - 13 H Left	4.0mm	13.4mm	184mm
	LCP668515	BW-015L	Locking Plate Anterolateral Distal Tibia 3.5mm - 15 H Left	4.0mm	13.4mm	210mm
80			LOCKING PLATE TIBIA DISTAL LATERAL 4.5/5.0mm, RIGHT/LEFT			
	LCP505102	TD-002R	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 2H Right	3.5mm	16.0mm	87mm
	LCP505104	TD-004R	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 4H Right	3.5mm	16.0mm	114mm
	LCP505106	TD-006R	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 6H Right	3.5mm	16.0mm	141mm
	LCP505108	TD-008R	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 8H Right	3.5mm	16.0mm	168mm
	LCP505110	TD-010R	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 10H Right	3.5mm	16.0mm	195mm
	LCP505202	TD-002L	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 2H Left	3.5mm	16.0mm	87mm
	LCP505204	TD-004L	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 4H Left	3.5mm	16.0mm	114mm
	LCP505206	TD-006L	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 6H Left	3.5mm	16.0mm	141mm
	LCP505208	TD-008L	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 8H Left	3.5mm	16.0mm	168mm
	LCP505210	TD-010L	Locking Plate Tibia Distal Lateral 4.5/5.0mm - 10H Left	3.5mm	16.0mm	195mm
81			LOCKING PLATE ANATOMICAL DISTAL TIBIA 3.5mm, RIGHT/LEFT			
	LCP677505	ADT-005R	Locking Plate Anatomical Distal Tibia 3.5mm - 5 H Right	4.0mm	13.2mm	116mm
	LCP677507	ADT-007R	Locking Plate Anatomical Distal Tibia 3.5mm - 7 H Right	4.0mm	13.2mm	142mm
	LCP677509	ADT-009R	Locking Plate Anatomical Distal Tibia 3.5mm - 9 H Right	4.0mm	13.2mm	168mm
	LCP677511	ADT-011R	Locking Plate Anatomical Distal Tibia 3.5mm - 11 H Right	4.0mm	13.2mm	194mm

	LCP677513	ADT-013R	Locking Plate Anatomical Distal Tibia 3.5mm - 13 H Right	4.0mm	13.2mm	220mm
	LCP677605	ADT-005L	Locking Plate Anatomical Distal Tibia 3.5mm - 5 H Left	4.0mm	13.2mm	116mm
	LCP677607	ADT-007L	Locking Plate Anatomical Distal Tibia 3.5mm - 7 H Left	4.0mm	13.2mm	142mm
	LCP677609	ADT-009L	Locking Plate Anatomical Distal Tibia 3.5mm - 9 H Left	4.0mm	13.2mm	168mm
	LCP677611	ADT-011L	Locking Plate Anatomical Distal Tibia 3.5mm - 11 H Left	4.0mm	13.2mm	194mm
	LCP677613	ADT-013L	Locking Plate Anatomical Distal Tibia 3.5mm - 13 H Left	4.0mm	13.2mm	220mm
82			LOCKING PLATE ANATOMICAL DISTAL TIBIA 4.5/5.0mm, RIGHT/LEFT			
	LCP678705	AQ-005R	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 5 H Right	4.5mm	14.5mm	142mm
	LCP678707	AQ-007R	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 7 H Right	4.5mm	14.5mm	168mm
	LCP678709	AQ-009R	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 9 H Right	4.5mm	14.5mm	194mm
	LCP678711	AQ-011R	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 11 H Right	4.5mm	14.5mm	220mm
	LCP678713	AQ-013R	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 13 H Right	4.5mm	14.5mm	246mm
	LCP678705	AQ-005L	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 5 H Right	4.5mm	14.5mm	142mm
	LCP678807	AQ-007L	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 7 H Left	4.5mm	14.5mm	168mm
	LCP678809	AQ-009L	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 9 H Left	4.5mm	14.5mm	194mm
	LCP678811	AQ-011L	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 11 H Left	4.5mm	14.5mm	220mm
	LCP678813	AQ-013L	Locking Plate Anatomical Distal Tibia 4.5/5.0mm - 13 H Left	4.5mm	14.5mm	246mm
83			LOCKING PLATE PILON 3.5mm			
	LCP669007	BS-007	Locking Plate Pilon 3.5mm - 7H	2.8mm	10.7mm	147mm
	LCP669009	BS-009	Locking Plate Pilon 3.5mm - 9H	2.8mm	10.7mm	173mm
84			LOCKING PLATE FIBULA DISTAL LATERAL 3.5mm			
	LCP235005	FD-005	Locking Plate Fibula Distal Lateral 3.5mm - 5H	2.8mm	10.0mm	64mm
	LCP235006	FD-006	Locking Plate Fibula Distal Lateral 3.5mm - 6H	2.8mm	10.0mm	82mm
	LCP235007	FD-007	Locking Plate Fibula Distal Lateral 3.5mm - 7H	2.8mm	10.0mm	100mm
	LCP235008	FD-008	Locking Plate Fibula Distal Lateral 3.5mm - 8H	2.8mm	10.0mm	118mm
	LCP235009	FD-009	Locking Plate Fibula Distal Lateral 3.5mm - 9H	2.8mm	10.0mm	136mm
	LCP235010	FD-010	Locking Plate Fibula Distal Lateral 3.5mm - 10H	2.8mm	10.0mm	154mm
	LCP235011	FD-011	Locking Plate Fibula Distal Lateral 3.5mm - 11H	2.8mm	10.0mm	172mm
	LCP235012	FD-012	Locking Plate Fibula Distal Lateral 3.5mm - 12H	2.8mm	10.0mm	190mm
	LCP235013	FD-013	Locking Plate Fibula Distal Lateral 3.5mm - 13H	2.8mm	10.0mm	208mm
	LCP235014	FD-014	Locking Plate Fibula Distal Lateral 3.5mm - 14H	2.8mm	10.0mm	226mm
85			LOCKING PLATE DISTAL FIBULA LATERAL 2.7/3.5mm, RIGHT/LEFT			
	LCP679603	LDFP-003R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 3H Right	3.5mm	12mm	73mm
	LCP679604	LDFP-004R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 4H Right	3.5mm	12mm	86mm
	LCP679605	LDFP-005R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 5H Right	3.5mm	12mm	99mm
	LCP679606	LDFP-006R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 6H Right	3.5mm	12mm	112mm

	LCP679607	LDFP-007R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 7H Right	3.5mm	12mm	125mm
	LCP679609	LDFP-009R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 9H Right	3.5mm	12mm	151mm
	LCP679611	LDFP-011R	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 11H Right	3.5mm	12mm	177mm
	LCP679503	LDFP-003L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 3H Left	3.5mm	12mm	73mm
	LCP679504	LDFP-004L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 4H Left	3.5mm	12mm	86mm
	LCP679505	LDFP-005L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 5H Left	3.5mm	12mm	99mm
	LCP679506	LDFP-006L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 6H Left	3.5mm	12mm	112mm
	LCP679507	LDFP-007L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 7H Left	3.5mm	12mm	125mm
	LCP679509	LDFP-009L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 9H Left	3.5mm	12mm	151mm
	LCP679511	LDFP-011L	Locking Plate Distal Fibula Lateral 2.7/3.5mm - 11H Left	3.5mm	12mm	177mm
86			LOCKING POSTEROLATERAL DISTAL FIBULA PLATE 2.7/3.5mm, RIGHT/LEFT			
	LCP679803	PDFP-003R	Locking Posterolateral Distal Fibula Plate - 3 H Right	3.5mm	12mm	77mm
	LCP679804	PDFP-004R	Locking Posterolateral Distal Fibula Plate - 4 H Right	3.5mm	12mm	90mm
	LCP679805	PDFP-005R	Locking Posterolateral Distal Fibula Plate - 5 H Right	3.5mm	12mm	103mm
	LCP679806	PDFP-006R	Locking Posterolateral Distal Fibula Plate - 6 H Right	3.5mm	12mm	116mm
	LCP679807	PDFP-007R	Locking Posterolateral Distal Fibula Plate - 7 H Right	3.5mm	12mm	129mm
	LCP679809	PDFP-009R	Locking Posterolateral Distal Fibula Plate - 9 H Right	3.5mm	12mm	155mm
	LCP679811	PDFP-011R	Locking Posterolateral Distal Fibula Plate - 11 H Right	3.5mm	12mm	181mm
	LCP679703	PDFP-003L	Locking Posterolateral Distal Fibula Plate - 3 H Left	3.5mm	12mm	77mm
	LCP679704	PDFP-004L	Locking Posterolateral Distal Fibula Plate - 4 H Left	3.5mm	12mm	90mm
	LCP679705	PDFP-005L	Locking Posterolateral Distal Fibula Plate - 5 H Left	3.5mm	12mm	103mm
	LCP679706	PDFP-006L	Locking Posterolateral Distal Fibula Plate - 6 H Left	3.5mm	12mm	116mm
	LCP679707	PDFP-007L	Locking Posterolateral Distal Fibula Plate - 7 H Left	3.5mm	12mm	129mm
	LCP679709	PDFP-009L	Locking Posterolateral Distal Fibula Plate - 9 H Left	3.5mm	12mm	155mm
	LCP679711	PDFP-011L	Locking Posterolateral Distal Fibula Plate - 11 H Left	3.5mm	12mm	181mm
87			LOCKING ANATOMICAL DISTAL TIBIA PLATE FOR FIBULA 3.5mm, RIGHT/LEFT			
	LCP679305	AT-005R	Locking Anatomical Distal Tibia Plate For Fibula - 5 H Right	3.0mm	13.5mm	116mm
	LCP679307	AT-007R	Locking Anatomical Distal Tibia Plate For Fibula - 7 H Right	3.0mm	13.5mm	142mm
	LCP679309	AT-009R	Locking Anatomical Distal Tibia Plate For Fibula - 9 H Right	3.0mm	13.5mm	168mm
	LCP679311	AT-011R	Locking Anatomical Distal Tibia Plate For Fibula - 11 H Right	3.0mm	13.5mm	194mm
	LCP679313	AT-013R	Locking Anatomical Distal Tibia Plate For Fibula - 13 H Right	3.0mm	13.5mm	220mm
	LCP679405	AT-005L	Locking Anatomical Distal Tibia Plate For Fibula - 5 H Left	3.0mm	13.5mm	116mm
	LCP679407	AT-007L	Locking Anatomical Distal Tibia Plate For Fibula - 7 H Left	3.0mm	13.5mm	142mm
	LCP679409	AT-009L	Locking Anatomical Distal Tibia Plate For Fibula - 9 H Left	3.0mm	13.5mm	168mm
	LCP679411	AT-011L	Locking Anatomical Distal Tibia Plate For Fibula - 11 H Left	3.0mm	13.5mm	194mm
	LCP679413	AT-013L	Locking Anatomical Distal Tibia Plate For Fibula - 13 H Left	3.0mm	13.5mm	220mm
88			LOCKING ANATOMICAL DISTAL TIBIA PLATE FOR FIBULA 4.5/5.0mm, RIGHT/LEFT			

	LCP678305	ATF-005R	Locking Anatomical Distal Tibia Plate For Fibula - 5 H Right	4.0mm	15.0mm	142mm
	LCP678307	ATF-007R	Locking Anatomical Distal Tibia Plate For Fibula - 7 H Right	4.0mm	15.0mm	168mm
	LCP678309	ATF-009R	Locking Anatomical Distal Tibia Plate For Fibula - 9 H Right	4.0mm	15.0mm	194mm
	LCP678311	ATF-011R	Locking Anatomical Distal Tibia Plate For Fibula - 11 H Right	4.0mm	15.0mm	220mm
	LCP678313	ATF-013R	Locking Anatomical Distal Tibia Plate For Fibula - 13 H Right	4.0mm	15.0mm	246mm
	LCP678405	ATF-005L	Locking Anatomical Distal Tibia Plate For Fibula - 5 H Left	4.0mm	15.0mm	142mm
	LCP678407	ATF-007L	Locking Anatomical Distal Tibia Plate For Fibula - 7 H Left	4.0mm	15.0mm	168mm
	LCP678409	ATF-009L	Locking Anatomical Distal Tibia Plate For Fibula - 9 H Left	4.0mm	15.0mm	194mm
	LCP678411	ATF-011L	Locking Anatomical Distal Tibia Plate For Fibula - 11 H Left	4.0mm	15.0mm	220mm
	LCP678413	ATF-013L	Locking Anatomical Distal Tibia Plate For Fibula - 13 H Left	4.0mm	15.0mm	246mm
89			LOCKING PLATE CALCANEAL 3.5mm, RIGHT/LEFT			
	LCP860113	AX-013R	Locking Plate Calcaneal 3.5mm - 13H Right	2.0mm		70mm
	LCP860115	AX-015R	Locking Plate Calcaneal 3.5mm - 15H Right	2.0mm		70mm
	LCP860013	AX-013L	Locking Plate Calcaneal 3.5mm - 13H Left	2.0mm		70mm
	LCP860015	AX-015L	Locking Plate Calcaneal 3.5mm - 15H Left	2.0mm		70mm
90			DHS/DCS Lag Screw			
	LCP600060	ZZ-060	DHS/DCS Lag Screw - 60mm	12.5mm		60mm
	LCP600065	ZZ-065	DHS/DCS Lag Screw - 65mm	12.5mm		65mm
	LCP600070	ZZ-070	DHS/DCS Lag Screw - 70mm	12.5mm		70mm
	LCP600075	ZZ-075	DHS/DCS Lag Screw - 75mm	12.5mm		75mm
	LCP600080	ZZ-080	DHS/DCS Lag Screw - 80mm	12.5mm		80mm
	LCP600085	ZZ-085	DHS/DCS Lag Screw - 85mm	12.5mm		85mm
	LCP600090	ZZ-090	DHS/DCS Lag Screw - 90mm	12.5mm		90mm
	LCP600095	ZZ-095	DHS/DCS Lag Screw - 95mm	12.5mm		95mm
	LCP600100	ZZ-100	DHS/DCS Lag Screw - 100mm	12.5mm		100mm
	LCP600105	ZZ-105	DHS/DCS Lag Screw - 105mm	12.5mm		105mm
	LCP600110	ZZ-110	DHS/DCS Lag Screw - 110mm	12.5mm		110mm
	LCP600115	ZZ-115	DHS/DCS Lag Screw - 115mm	12.5mm		115mm
	LCP600120	ZZ-120	DHS/DCS Lag Screw - 120mm	12.5mm		120mm
			COMPRESSION SCREW			
	LCP600000	ZZ-000	Compression Screw - 36mm	3.8mm		36mm
91			CORTICAL SCREW, Dia - 2.4mm - SELF TAPPING			
	701.06ST-Ti	701.06ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 06mm	2.4mm		06mm
	701.08ST-Ti	701.08ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 08mm	2.4mm		08mm
	701.10ST-Ti	701.10ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 10mm	2.4mm		10mm
	701.12ST-Ti	701.12ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 12mm	2.4mm		12mm

	701.14ST-Ti	701.14ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 14mm	2.4mm		14mm
	701.16ST-Ti	701.16ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 16mm	2.4mm		16mm
	701.18ST-Ti	701.18ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 18mm	2.4mm		18mm
	701.20ST-Ti	701.20ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 20mm	2.4mm		20mm
	701.22ST-Ti	701.22ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 22mm	2.4mm		22mm
	701.24ST-Ti	701.24ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 24mm	2.4mm		24mm
	701.26ST-Ti	701.26ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 26mm	2.4mm		26mm
	701.28ST-Ti	701.28ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 28mm	2.4mm		28mm
	701.30ST-Ti	701.30ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 30mm	2.4mm		30mm
	701.32ST-Ti	701.32ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 32mm	2.4mm		32mm
	701.34ST-Ti	701.34ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 34mm	2.4mm		34mm
	701.36ST-Ti	701.36ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 36mm	2.4mm		36mm
	701.38ST-Ti	701.38ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 38mm	2.4mm		38mm
	701.40ST-Ti	701.40ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 40mm	2.4mm		40mm
	701.42ST-Ti	701.42ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 42mm	2.4mm		42mm
	701.44ST-Ti	701.44ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 44mm	2.4mm		44mm
	701.46ST-Ti	701.46ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 46mm	2.4mm		46mm
	701.48ST-Ti	701.48ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 48mm	2.4mm		48mm
	701.50ST-Ti	701.50ST	Cortical Screw Dia - 2.4mm - Self Tapping Length - 50mm	2.4mm		50mm
92			CORTICAL SCREW, Dia - 2.7mm - SELF TAPPING			
	603.06ST-Ti	603.06ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 06mm	2.7mm		06mm
	603.08ST-Ti	603.08ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 08mm	2.7mm		08mm
	603.10ST-Ti	603.10ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 10mm	2.7mm		10mm
	603.12ST-Ti	603.12ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 12mm	2.7mm		12mm
	603.14ST-Ti	603.14ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 14mm	2.7mm		14mm
	603.16ST-Ti	603.16ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 16mm	2.7mm		16mm
	603.18ST-Ti	603.18ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 18mm	2.7mm		18mm
	603.20ST-Ti	603.20ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 20mm	2.7mm		20mm
	603.22ST-Ti	603.22ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 22mm	2.7mm		22mm
	603.24ST-Ti	603.24ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 24mm	2.7mm		24mm
	603.26ST-Ti	603.26ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 26mm	2.7mm		26mm
	603.28ST-Ti	603.28ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 28mm	2.7mm		28mm
	603.30ST-Ti	603.30ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 30mm	2.7mm		30mm
	603.32ST-Ti	603.32ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 32mm	2.7mm		32mm
	603.34ST-Ti	603.34ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 34mm	2.7mm		34mm
	603.36ST-Ti	603.36ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 36mm	2.7mm		36mm
	603.38ST-Ti	603.38ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 38mm	2.7mm		38mm
	603.40ST-Ti	603.40ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 40mm	2.7mm		40mm
	603.42ST-Ti	603.42ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 42mm	2.7mm		42mm
	603.44ST-Ti	603.44ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 44mm	2.7mm		44mm
	603.46ST-Ti	603.46ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 46mm	2.7mm		46mm
	603.48ST-Ti	603.48ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 48mm	2.7mm		48mm

	603.50ST-Ti	603.50ST	Cortical Screw Dia - 2.7mm - Self Tapping Length - 50mm	2.7mm		50mm
93			LOCKING HEAD SCREW, Dia - 2.7mm - SELF TAPPING			
	LCP580006	DD-006	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 6mm	2.7mm		06mm
	LCP580008	DD-008	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 8mm	2.7mm		08mm
	LCP580010	DD-010	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 10mm	2.7mm		10mm
	LCP580012	DD-012	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 12mm	2.7mm		12mm
	LCP580014	DD-014	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 14mm	2.7mm		14mm
	LCP580016	DD-016	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 16mm	2.7mm		16mm
	LCP580018	DD-018	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 18mm	2.7mm		18mm
	LCP580020	DD-020	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 20mm	2.7mm		20mm
	LCP580022	DD-022	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 22mm	2.7mm		22mm
	LCP580024	DD-024	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 24mm	2.7mm		24mm
	LCP580026	DD-026	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 26mm	2.7mm		26mm
	LCP580028	DD-028	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 28mm	2.7mm		28mm
	LCP580030	DD-030	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 30mm	2.7mm		30mm
	LCP580032	DD-032	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 32mm	2.7mm		32mm
	LCP580034	DD-034	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 34mm	2.7mm		34mm
	LCP580036	DD-036	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 36mm	2.7mm		36mm
	LCP580038	DD-038	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 38mm	2.7mm		38mm
	LCP580040	DD-040	Locking Head Screw, Dia - 2.7mm - Self Tapping Length - 40mm	2.7mm		40mm
94			CORTICAL SCREW, Dia - 3.5mm - SELF TAPPING			
	604.14ST-Ti	604.14ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 14mm	3.5mm		14mm
	604.16ST-Ti	604.16ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 16mm	3.5mm		16mm
	604.18ST-Ti	604.18ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 18mm	3.5mm		18mm
	604.20ST-Ti	604.20ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 20mm	3.5mm		20mm
	604.22ST-Ti	604.22ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 22mm	3.5mm		22mm
	604.24ST-Ti	604.24ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 24mm	3.5mm		24mm
	604.26ST-Ti	604.26ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 26mm	3.5mm		26mm
	604.28ST-Ti	604.28ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 28mm	3.5mm		28mm
	604.30ST-Ti	604.30ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 30mm	3.5mm		30mm
	604.32ST-Ti	604.32ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 32mm	3.5mm		32mm
	604.34ST-Ti	604.34ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 34mm	3.5mm		34mm
	604.36ST-Ti	604.36ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 36mm	3.5mm		36mm
	604.38ST-Ti	604.38ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 38mm	3.5mm		38mm
	604.40ST-Ti	604.40ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 40mm	3.5mm		40mm
	604.45ST-Ti	604.45ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 45mm	3.5mm		45mm
	604.50ST-Ti	604.50ST	Cortical Screw, Dia - 3.5mm - Self Tapping Length - 50mm	3.5mm		50mm
95			CORTICAL SCREW, Dia - 4.5mm - SELF TAPPING			
	605.22ST-Ti	605.22ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 22mm	4.5mm		22mm
	605.24ST-Ti	605.24ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 24mm	4.5mm		24mm

	605.26ST-Ti	605.26ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 26mm	4.5mm		26mm
	605.28ST-Ti	605.28ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 28mm	4.5mm		28mm
	605.30ST-Ti	605.30ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 30mm	4.5mm		30mm
	605.32ST-Ti	605.32ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 32mm	4.5mm		32mm
	605.34ST-Ti	605.34ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 34mm	4.5mm		34mm
	605.36ST-Ti	605.36ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 36mm	4.5mm		36mm
	605.38ST-Ti	605.38ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 38mm	4.5mm		38mm
	605.40ST-Ti	605.40ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 40mm	4.5mm		40mm
	605.42ST-Ti	605.42ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 42mm	4.5mm		42mm
	605.44ST-Ti	605.44ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 44mm	4.5mm		44mm
	605.46ST-Ti	605.46ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 46mm	4.5mm		46mm
	605.48ST-Ti	605.48ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 48mm	4.5mm		48mm
	605.50ST-Ti	605.50ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 50mm	4.5mm		50mm
	605.52ST-Ti	605.52ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 52mm	4.5mm		52mm
	605.54ST-Ti	605.54ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 54mm	4.5mm		54mm
	605.56ST-Ti	605.56ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 56mm	4.5mm		56mm
	605.58ST-Ti	605.58ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 58mm	4.5mm		58mm
	605.60ST-Ti	605.60ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 60mm	4.5mm		60mm
	605.62ST-Ti	605.62ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 62mm	4.5mm		62mm
	605.64ST-Ti	605.64ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 64mm	4.5mm		64mm
	605.66ST-Ti	605.66ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 66mm	4.5mm		66mm
	605.68ST-Ti	605.68ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 68mm	4.5mm		68mm
	605.70ST-Ti	605.70ST	Cortical Screw, Dia - 4.5mm - Self Tapping Length - 70mm	4.5mm		70mm
96			LOCKING HEAD SCREW, Dia - 3.5mm - SELF TAPPING			
	LCP451010	OO-010	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 10mm	3.5mm		10mm
	LCP451012	OO-012	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 12mm	3.5mm		12mm
	LCP451014	OO-014	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 14mm	3.5mm		14mm
	LCP451016	OO-016	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 16mm	3.5mm		16mm
	LCP451018	OO-018	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 18mm	3.5mm		18mm
	LCP451020	OO-020	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 20mm	3.5mm		20mm
	LCP451022	OO-022	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 22mm	3.5mm		22mm
	LCP451024	OO-024	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 24mm	3.5mm		24mm
	LCP451026	OO-026	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 26mm	3.5mm		26mm
	LCP451028	OO-028	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 28mm	3.5mm		28mm
	LCP451030	OO-030	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 30mm	3.5mm		30mm
	LCP451032	OO-032	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 32mm	3.5mm		32mm
	LCP451034	OO-034	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 34mm	3.5mm		34mm
	LCP451036	OO-036	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 36mm	3.5mm		36mm
	LCP451038	OO-038	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 38mm	3.5mm		38mm
	LCP451040	OO-040	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 40mm	3.5mm		40mm
	LCP451045	OO-045	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 45mm	3.5mm		45mm
	LCP451050	OO-050	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 50mm	3.5mm		50mm

	LCP451055	OO-055	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 55mm	3.5mm	55mm
	LCP451060	OO-060	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 60mm	3.5mm	60mm
	LCP451065	OO-065	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 65mm	3.5mm	65mm
	LCP451070	OO-070	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 70mm	3.5mm	70mm
	LCP451075	OO-075	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 75mm	3.5mm	75mm
	LCP451080	OO-080	Locking Head Screw, Dia - 3.5mm - Self Tapping Length - 80mm	3.5mm	80mm
97			LOCKING HEAD SCREW, Dia - 5.0mm - SELF TAPPING		
	LCP552014	PP-014	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 14mm	5.0mm	14mm
	LCP552016	PP-016	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 16mm	5.0mm	16mm
	LCP552018	PP-018	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 18mm	5.0mm	18mm
	LCP552020	PP-020	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 20mm	5.0mm	20mm
	LCP552022	PP-022	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 22mm	5.0mm	22mm
	LCP552024	PP-024	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 24mm	5.0mm	24mm
	LCP552026	PP-026	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 26mm	5.0mm	26mm
	LCP552028	PP-028	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 28mm	5.0mm	28mm
	LCP552030	PP-030	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 30mm	5.0mm	30mm
	LCP552032	PP-032	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 32mm	5.0mm	32mm
	LCP552034	PP-034	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 34mm	5.0mm	34mm
	LCP552036	PP-036	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 36mm	5.0mm	36mm
	LCP552038	PP-038	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 38mm	5.0mm	38mm
	LCP552040	PP-040	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 40mm	5.0mm	40mm
	LCP552042	PP-042	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 42mm	5.0mm	42mm
	LCP552044	PP-044	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 44mm	5.0mm	44mm
	LCP552046	PP-046	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 46mm	5.0mm	46mm
	LCP552048	PP-048	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 48mm	5.0mm	48mm
	LCP552050	PP-050	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 50mm	5.0mm	50mm
	LCP552055	PP-055	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 55mm	5.0mm	55mm
	LCP552060	PP-060	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 60mm	5.0mm	60mm
	LCP552065	PP-065	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 65mm	5.0mm	65mm
	LCP552070	PP-070	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 70mm	5.0mm	70mm
	LCP552075	PP-075	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 75mm	5.0mm	75mm
	LCP552080	PP-080	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 80mm	5.0mm	80mm
	LCP552085	PP-085	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 85mm	5.0mm	85mm
	LCP552090	PP-090	Locking Head Screw, Dia - 5.0mm - Self Tapping Length - 90mm	5.0mm	90mm
98			LOCKING HEAD CANCELLOUS SCREW Dia - 5.0mm "16mm THREAD"		
	LCP501630	5016-30	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 30mm	5.0mm	30mm
	LCP501635	5016-35	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 35mm	5.0mm	35mm
	LCP501640	5016-40	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 40mm	5.0mm	40mm
	LCP501645	5016-45	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 45mm	5.0mm	45mm
	LCP501650	5016-50	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 50mm	5.0mm	50mm
	LCP501655	5016-55	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 55mm	5.0mm	55mm

	LCP501660	5016-60	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 60mm	5.0mm	60mm
	LCP501665	5016-65	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 65mm	5.0mm	65mm
	LCP501670	5016-70	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 70mm	5.0mm	70mm
	LCP501675	5016-75	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 75mm	5.0mm	75mm
	LCP501680	5016-80	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 80mm	5.0mm	80mm
	LCP501685	5016-85	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 85mm	5.0mm	85mm
	LCP501690	5016-90	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 90mm	5.0mm	90mm
	LCP501695	5016-95	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 95mm	5.0mm	95mm
	LCP5016100	5016-100	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 100mm	5.0mm	100mm
	LCP5016105	5016-105	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 105mm	5.0mm	105mm
	LCP5016110	5016-110	Locking Head Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 110mm	5.0mm	110mm
99			LOCKING HEAD CANCELLOUS SCREW Dia - 5.0mm "32mm THREAD"		
	LCP503240	5032-40	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 40mm	5.0mm	40mm
	LCP503245	5032-45	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 45mm	5.0mm	45mm
	LCP503250	5032-50	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 50mm	5.0mm	50mm
	LCP503255	5032-55	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 55mm	5.0mm	55mm
	LCP503260	5032-60	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 60mm	5.0mm	60mm
	LCP503265	5032-65	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 65mm	5.0mm	65mm
	LCP503270	5032-70	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 70mm	5.0mm	70mm
	LCP503275	5032-75	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 75mm	5.0mm	75mm
	LCP503280	5032-80	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 80mm	5.0mm	80mm
	LCP503285	5032-85	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 85mm	5.0mm	85mm
	LCP503290	5032-90	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 90mm	5.0mm	90mm
	LCP503295	5032-95	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 95mm	5.0mm	95mm
	LCP5032100	5032-100	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 100mm	5.0mm	100mm
	LCP5032105	5032-105	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 105mm	5.0mm	105mm
	LCP5032110	5032-110	Locking Head Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 110mm	5.0mm	110mm
100			LOCKING HEAD CANCELLOUS SCREW Dia - 5.0mm "FULL THREAD"		
	LCP500025	5000-25	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 25mm	5.0mm	25mm
	LCP500030	5000-30	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 30mm	5.0mm	30mm
	LCP500035	5000-35	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 35mm	5.0mm	35mm
	LCP500040	5000-40	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 40mm	5.0mm	40mm
	LCP500045	5000-45	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 45mm	5.0mm	45mm
	LCP500050	5000-50	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 50mm	5.0mm	50mm
	LCP500055	5000-55	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 55mm	5.0mm	55mm
	LCP500060	5000-60	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 60mm	5.0mm	60mm
	LCP500065	5000-65	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 65mm	5.0mm	65mm
	LCP500070	5000-70	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 70mm	5.0mm	70mm
	LCP500075	5000-75	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 75mm	5.0mm	75mm
	LCP500080	5000-80	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 80mm	5.0mm	80mm
	LCP500085	5000-85	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 85mm	5.0mm	85mm

	LCP500090	5000-90	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 90mm	5.0mm		90mm
	LCP500095	5000-95	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 95mm	5.0mm		95mm
	LCP5000100	5000-100	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 100mm	5.0mm		100mm
	LCP5000105	5000-105	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 105mm	5.0mm		105mm
	LCP5000110	5000-110	Locking Head Cancellous Screw Dia - 5.0mm "Full Thread" Length - 110mm	5.0mm		110mm
101			LOCKING HEAD CANCELLOUS SCREW Dia - 6.5mm "16mm THREAD"			
	LCP651630	6516-30	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 30mm	6.5mm		30mm
	LCP651635	6516-35	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 35mm	6.5mm		35mm
	LCP651640	6516-40	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 40mm	6.5mm		40mm
	LCP651645	6516-45	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 45mm	6.5mm		45mm
	LCP651650	6516-50	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 50mm	6.5mm		50mm
	LCP651655	6516-55	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 55mm	6.5mm		55mm
	LCP651660	6516-60	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 60mm	6.5mm		60mm
	LCP651665	6516-65	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 65mm	6.5mm		65mm
	LCP651670	6516-70	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 70mm	6.5mm		70mm
	LCP651675	6516-75	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 75mm	6.5mm		75mm
	LCP651680	6516-80	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 80mm	6.5mm		80mm
	LCP651685	6516-85	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 85mm	6.5mm		85mm
	LCP651690	6516-90	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 90mm	6.5mm		90mm
	LCP651695	6516-95	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 95mm	6.5mm		95mm
	LCP6516100	6516-100	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 100mm	6.5mm		100mm
	LCP6516105	6516-105	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 105mm	6.5mm		105mm
	LCP6516110	6516-110	Locking Head Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 110mm	6.5mm		110mm
102			LOCKING HEAD CANCELLOUS SCREW Dia - 6.5mm "32mm THREAD"			
	LCP653240	6532-40	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 40mm	6.5mm		40mm
	LCP653245	6532-45	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 45mm	6.5mm		45mm
	LCP653250	6532-50	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 50mm	6.5mm		50mm
	LCP653255	6532-55	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 55mm	6.5mm		55mm
	LCP653260	6532-60	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 60mm	6.5mm		60mm
	LCP653265	6532-65	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 65mm	6.5mm		65mm
	LCP653270	6532-70	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 70mm	6.5mm		70mm
	LCP653275	6532-75	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 75mm	6.5mm		75mm
	LCP653280	6532-80	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 80mm	6.5mm		80mm
	LCP653285	6532-85	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 85mm	6.5mm		85mm
	LCP653290	6532-90	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 90mm	6.5mm		90mm
	LCP653295	6532-95	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 95mm	6.5mm		95mm
	LCP6532100	6532-100	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 100mm	6.5mm		100mm
	LCP6532105	6532-105	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 105mm	6.5mm		105mm
	LCP6532110	6532-110	Locking Head Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 110mm	6.5mm		110mm
103			LOCKING HEAD CANCELLOUS SCREW Dia - 6.5mm "FULL THREAD"			

	LCP650025	6500-25	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 25mm	6.5mm	25mm
	LCP650030	6500-30	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 30mm	6.5mm	30mm
	LCP650035	6500-35	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 35mm	6.5mm	35mm
	LCP650040	6500-40	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 40mm	6.5mm	40mm
	LCP650045	6500-45	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 45mm	6.5mm	45mm
	LCP650050	6500-50	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 50mm	6.5mm	50mm
	LCP650055	6500-55	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 55mm	6.5mm	55mm
	LCP650060	6500-60	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 60mm	6.5mm	60mm
	LCP650065	6500-65	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 65mm	6.5mm	65mm
	LCP650070	6500-70	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 70mm	6.5mm	70mm
	LCP650075	6500-75	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 75mm	6.5mm	75mm
	LCP650080	6500-80	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 80mm	6.5mm	80mm
	LCP650085	6500-85	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 85mm	6.5mm	85mm
	LCP650090	6500-90	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 90mm	6.5mm	90mm
	LCP650095	6500-95	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 95mm	6.5mm	95mm
	LCP650100	6500-100	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 100mm	6.5mm	100mm
	LCP650105	6500-105	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 105mm	6.5mm	105mm
	LCP650110	6500-110	Locking Head Cancellous Screw Dia - 6.5mm "Full Thread" Length - 110mm	6.5mm	110mm
104			Locking Head Cancellous Screw 3.5 mm "Full Thread"		
	LCP-35F-014	35F-014	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 14mm	3.5mm	14mm
	LCP-35F-016	35F-016	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 16mm	3.5mm	16mm
	LCP-35F-018	35F-018	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 18mm	3.5mm	18mm
	LCP-35F-020	35F-020	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 20mm	3.5mm	20mm
	LCP-35F-022	35F-022	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 22mm	3.5mm	22mm
	LCP-35F-024	35F-024	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 24mm	3.5mm	24mm
	LCP-35F-026	35F-026	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 26mm	3.5mm	26mm
	LCP-35F-028	35F-028	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 28mm	3.5mm	28mm
	LCP-35F-030	35F-030	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 30mm	3.5mm	30mm
	LCP-35F-035	35F-035	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 35mm	3.5mm	35mm
	LCP-35F-040	35F-040	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 40mm	3.5mm	40mm
	LCP-35F-045	35F-045	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 45mm	3.5mm	45mm
	LCP-35F-050	35F-050	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 50mm	3.5mm	50mm
	LCP-35F-055	35F-055	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 55mm	3.5mm	55mm
	LCP-35F-060	35F-060	Locking Head Cancellous Screw 3.5mm "Full Thread" Length - 60mm	3.5mm	60mm
105			Locking Head Cancellous Screw 3.5 mm "Short Thread"		
	LCP-35S-014	35S-014	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 14mm	3.5mm	14mm
	LCP-35S-016	35S-016	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 16mm	3.5mm	16mm
	LCP-35S-018	35S-018	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 18mm	3.5mm	18mm
	LCP-35S-020	35S-020	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 20mm	3.5mm	20mm
	LCP-35S-022	35S-022	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 22mm	3.5mm	22mm
	LCP-35S-024	35S-024	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 24mm	3.5mm	24mm

	LCP-35S-026	35S-026	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 26mm	3.5mm		26mm
	LCP-35S-028	35S-028	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 28mm	3.5mm		28mm
	LCP-35S-030	35S-030	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 30mm	3.5mm		30mm
	LCP-35S-035	35S-035	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 35mm	3.5mm		35mm
	LCP-35S-040	35S-040	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 40mm	3.5mm		40mm
	LCP-35S-045	35S-045	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 45mm	3.5mm		45mm
	LCP-35S-050	35S-050	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 50mm	3.5mm		50mm
	LCP-35S-055	35S-055	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 55mm	3.5mm		55mm
	LCP-35S-060	35S-060	Locking Head Cancellous Screw 3.5mm "Short Thread" Length - 60mm	3.5mm		60mm
106			LOCKING HEAD CANNULATED CANCELLOUS SCREW Dia - 5.0mm "16mm THREAD"			
	LCP511630	5116-30	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 30mm	5.0mm		30mm
	LCP511635	5116-35	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 35mm	5.0mm		35mm
	LCP511640	5116-40	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 40mm	5.0mm		40mm
	LCP511645	5116-45	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 45mm	5.0mm		45mm
	LCP511650	5116-50	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 50mm	5.0mm		50mm
	LCP511655	5116-55	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 55mm	5.0mm		55mm
	LCP511660	5116-60	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 60mm	5.0mm		60mm
	LCP511665	5116-65	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 65mm	5.0mm		65mm
	LCP511670	5116-70	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 70mm	5.0mm		70mm
	LCP511675	5116-75	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 75mm	5.0mm		75mm
	LCP511680	5116-80	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 80mm	5.0mm		80mm
	LCP511685	5116-85	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 85mm	5.0mm		85mm
	LCP511690	5116-90	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 90mm	5.0mm		90mm
	LCP511695	5116-95	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 95mm	5.0mm		95mm
	LCP5116100	5116-100	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 100mm	5.0mm		100mm
	LCP5116105	5116-105	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 105mm	5.0mm		105mm
	LCP5116110	5116-110	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "16mm Thread" Length - 110mm	5.0mm		110mm
107			LOCKING HEAD CANNULATED CANCELLOUS SCREW Dia - 5.0mm "32mm THREAD"			
	LCP513240	5132-40	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 40mm	5.0mm		40mm
	LCP513245	5132-45	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 45mm	5.0mm		45mm
	LCP513250	5132-50	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 50mm	5.0mm		50mm
	LCP513255	5132-55	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 55mm	5.0mm		55mm
	LCP513260	5132-60	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 60mm	5.0mm		60mm
	LCP513265	5132-65	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 65mm	5.0mm		65mm
	LCP513270	5132-70	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 70mm	5.0mm		70mm
	LCP513275	5132-75	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 75mm	5.0mm		75mm
	LCP513280	5132-80	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 80mm	5.0mm		80mm
	LCP513285	5132-85	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 85mm	5.0mm		85mm
	LCP513290	5132-90	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 90mm	5.0mm		90mm
	LCP513295	5132-95	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 95mm	5.0mm		95mm
	LCP5132100	5132-100	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 100mm	5.0mm		100mm

	LCP5132105	5132-105	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 105mm	5.0mm		105mm
	LCP5132110	5132-110	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "32mm Thread" Length - 110mm	5.0mm		110mm
108			LOCKING HEAD CANNULATED CANCELLOUS SCREW Dia - 5.0mm "FULL THREAD"			
	LCP510025	5100-25	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 25mm	5.0mm		25mm
	LCP510030	5100-30	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 30mm	5.0mm		30mm
	LCP510035	5100-35	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 35mm	5.0mm		35mm
	LCP510040	5100-40	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 40mm	5.0mm		40mm
	LCP510045	5100-45	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 45mm	5.0mm		45mm
	LCP510050	5100-50	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 50mm	5.0mm		50mm
	LCP510055	5100-55	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 55mm	5.0mm		55mm
	LCP510060	5100-60	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 60mm	5.0mm		60mm
	LCP510065	5100-65	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 65mm	5.0mm		65mm
	LCP510070	5100-70	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 70mm	5.0mm		70mm
	LCP510075	5100-75	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 75mm	5.0mm		75mm
	LCP510080	5100-80	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 80mm	5.0mm		80mm
	LCP510085	5100-85	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 85mm	5.0mm		85mm
	LCP510090	5100-90	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 90mm	5.0mm		90mm
	LCP510095	5100-95	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 95mm	5.0mm		95mm
	LCP5100100	5100-100	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 100mm	5.0mm		100mm
	LCP5100105	5100-105	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 105mm	5.0mm		105mm
	LCP5100110	5100-110	Locking Head Cannulated Cancellous Screw Dia - 5.0mm "Full Thread" Length - 110mm	5.0mm		110mm
109			Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread"			
	LCCS-35F-014	LC683.14	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 14mm Length	3.5mm		14mm
	LCCS-35F-016	LC683.16	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 16mm Length	3.5mm		16mm
	LCCS-35F-018	LC683.18	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 18mm Length	3.5mm		18mm
	LCCS-35F-020	LC683.20	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 20mm Length	3.5mm		20mm
	LCCS-35F-022	LC683.22	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 22mm Length	3.5mm		22mm
	LCCS-35F-024	LC683.24	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 24mm Length	3.5mm		24mm
	LCCS-35F-026	LC683.26	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 26mm Length	3.5mm		26mm
	LCCS-35F-028	LC683.28	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 28mm Length	3.5mm		28mm
	LCCS-35F-030	LC683.30	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 30mm Length	3.5mm		30mm
	LCCS-35F-035	LC683.35	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 35mm Length	3.5mm		35mm
	LCCS-35F-040	LC683.40	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 40mm Length	3.5mm		40mm
	LCCS-35F-045	LC683.45	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 45mm Length	3.5mm		45mm
	LCCS-35F-050	LC683.50	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 50mm Length	3.5mm		50mm
	LCCS-35F-055	LC683.55	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 55mm Length	3.5mm		55mm
	LCCS-35F-060	LC683.60	Locking Head Cannulated Cancellous Screw 3.5 mm "Full Thread" - 60mm Length	3.5mm		60mm
110			Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread"			
	LCCS-35S-014	LC684.14	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 14mm Length	3.5mm		14mm
	LCCS-35S-016	LC684.16	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 16mm Length	3.5mm		16mm

	LCCS-35S-018	LC684.18	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 18mm Length	3.5mm	18mm
	LCCS-35S-020	LC684.20	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 20mm Length	3.5mm	20mm
	LCCS-35S-022	LC684.22	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 22mm Length	3.5mm	22mm
	LCCS-35S-024	LC684.24	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 24mm Length	3.5mm	24mm
	LCCS-35S-026	LC684.26	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 26mm Length	3.5mm	26mm
	LCCS-35S-028	LC684.28	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 28mm Length	3.5mm	28mm
	LCCS-35S-030	LC684.30	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 30mm Length	3.5mm	30mm
	LCCS-35S-035	LC684.35	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 35mm Length	3.5mm	35mm
	LCCS-35S-040	LC684.40	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 40mm Length	3.5mm	40mm
	LCCS-35S-045	LC684.45	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 45mm Length	3.5mm	45mm
	LCCS-35S-050	LC684.50	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 50mm Length	3.5mm	50mm
	LCCS-35S-055	LC684.55	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 55mm Length	3.5mm	55mm
	LCCS-35S-060	LC684.60	Locking Head Cannulated Cancellous Screw 3.5 mm "Short Thread" - 60mm Length	3.5mm	60mm
111			LOCKING HEAD CANNULATED CANCELLOUS SCREW Dia - 6.5mm "16mm THREAD"		
	LCP661630	6116-30	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 30mm	6.5mm	30mm
	LCP661635	6116-35	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 35mm	6.5mm	35mm
	LCP661640	6116-40	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 40mm	6.5mm	40mm
	LCP661645	6116-45	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 45mm	6.5mm	45mm
	LCP661650	6116-50	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 50mm	6.5mm	50mm
	LCP661655	6116-55	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 55mm	6.5mm	55mm
	LCP661660	6116-60	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 60mm	6.5mm	60mm
	LCP661665	6116-65	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 65mm	6.5mm	65mm
	LCP661670	6116-70	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 70mm	6.5mm	70mm
	LCP661675	6116-75	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 75mm	6.5mm	75mm
	LCP661680	6116-80	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 80mm	6.5mm	80mm
	LCP661685	6116-85	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 85mm	6.5mm	85mm
	LCP661690	6116-90	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 90mm	6.5mm	90mm
	LCP661695	6116-95	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 95mm	6.5mm	95mm
	LCP6616100	6116-100	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 100mm	6.5mm	100mm
	LCP6616105	6116-105	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 105mm	6.5mm	105mm
	LCP6616110	6116-110	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "16mm Thread" Length - 110mm	6.5mm	110mm
112			LOCKING HEAD CANNULATED CANCELLOUS SCREW Dia - 6.5mm "32mm THREAD"		
	LCP663240	6132-40	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 40mm	6.5mm	40mm
	LCP663245	6132-45	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 45mm	6.5mm	45mm
	LCP663250	6132-50	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 50mm	6.5mm	50mm
	LCP663255	6132-55	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 55mm	6.5mm	55mm
	LCP663260	6132-60	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 60mm	6.5mm	60mm
	LCP663265	6132-65	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 65mm	6.5mm	65mm
	LCP663270	6132-70	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 70mm	6.5mm	70mm
	LCP663275	6132-75	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 75mm	6.5mm	75mm
	LCP663280	6132-80	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 80mm	6.5mm	80mm

	LCP663285	6132-85	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 85mm	6.5mm	85mm
	LCP663290	6132-90	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 90mm	6.5mm	90mm
	LCP663295	6132-95	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 95mm	6.5mm	95mm
	LCP6632100	6132-100	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 100mm	6.5mm	100mm
	LCP6632105	6132-105	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 105mm	6.5mm	105mm
	LCP6632110	6132-110	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "32mm Thread" Length - 110mm	6.5mm	110mm
113			LOCKING HEAD CANNULATED CANCELLOUS SCREW Dia - 6.5mm "FULL THREAD"		
	LCP660025	6100-25	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 25mm	6.5mm	25mm
	LCP660030	6100-30	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 30mm	6.5mm	30mm
	LCP660035	6100-35	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 35mm	6.5mm	35mm
	LCP660040	6100-40	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 40mm	6.5mm	40mm
	LCP660045	6100-45	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 45mm	6.5mm	45mm
	LCP660050	6100-50	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 50mm	6.5mm	50mm
	LCP660055	6100-55	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 55mm	6.5mm	55mm
	LCP660060	6100-60	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 60mm	6.5mm	60mm
	LCP660065	6100-65	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 65mm	6.5mm	65mm
	LCP660070	6100-70	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 70mm	6.5mm	70mm
	LCP660075	6100-75	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 75mm	6.5mm	75mm
	LCP660080	6100-80	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 80mm	6.5mm	80mm
	LCP660085	6100-85	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 85mm	6.5mm	85mm
	LCP660090	6100-90	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 90mm	6.5mm	90mm
	LCP660095	6100-95	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 95mm	6.5mm	95mm
	LCP660100	6100-100	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 100mm	6.5mm	100mm
	LCP6600105	6100-105	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 105mm	6.5mm	105mm
	LCP6600110	6100-110	Locking Head Cannulated Cancellous Screw Dia - 6.5mm "Full Thread" Length - 110mm	6.5mm	110mm
114			LOCKING HEAD SCREW, Dia - 2.4mm - SELF TAPPING		
	LCP570006	FF-006	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 6mm	2.4mm	06mm
	LCP570008	FF-008	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 8mm	2.4mm	08mm
	LCP570010	FF-010	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 10mm	2.4mm	10mm
	LCP570012	FF-012	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 12mm	2.4mm	12mm
	LCP570014	FF-014	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 14mm	2.4mm	14mm
	LCP570016	FF-016	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 16mm	2.4mm	16mm
	LCP570018	FF-018	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 18mm	2.4mm	18mm
	LCP570020	FF-020	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 20mm	2.4mm	20mm
	LCP570022	FF-022	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 22mm	2.4mm	22mm
	LCP570024	FF-024	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 24mm	2.4mm	24mm
	LCP570026	FF-026	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 26mm	2.4mm	26mm
	LCP570028	FF-028	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 28mm	2.4mm	28mm
	LCP570030	FF-030	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 30mm	2.4mm	30mm
	LCP570032	FF-032	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 32mm	2.4mm	32mm
	LCP570034	FF-034	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 34mm	2.4mm	34mm

	LCP570036	FF-036	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 36mm	2.4mm		36mm
	LCP570038	FF-038	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 38mm	2.4mm		38mm
	LCP570040	FF-040	Locking Head Screw, Dia - 2.4mm - Self Tapping Length - 40mm	2.4mm		40mm